

# HPLC

## Column Protection Guide

Version 0610

### **Includes:**

- Mobile Phase Limitations
- Column Storage Tips
- Column Protection Devices

# COLUMN PROTECTION GUIDE

## TABLE OF CONTENTS

### **INTRODUCTION TO COLUMN CARE AND USAGE 1-4**

Selecting the Right Tubing and Fittings  
Column Installation

### **PART I - SILICA-BASED & TWIN™ TECHNOLOGY COLUMNS 4-7**

Running Parameters  
Mobile Phase Considerations  
Stationary Phase Considerations  
Backpressure and Flow Rates

Scaling Up/Scaling Down  
Column Storage  
Column Cleaning Procedures

### **PART II - SFC (SUPERCRITICAL FLUID CHROMATOGRAPHY) 7**

Running Parameters  
Equilibrating Column  
Mobile Phase Considerations

Cleaning Procedure  
Column Storage

### **PART III - AXIA PREPARATIVE COLUMNS 8**

Running Parameters  
Mobile Phase Considerations

Cleaning Procedure  
Column Storage

### **PART IV - LUX CHIRAL COLUMNS 8-9**

Running Parameters  
Mobile Phase Considerations

Extending Lifetime and Reconditioning  
Column Storage

### **PART V - CHIREX CHIRAL COLUMNS 10**

Running Parameters  
Mobile Phase Considerations

### **PART VI - BIOSEP-SEC-S COLUMNS 10**

Running Parameters  
Mobile Phase Considerations  
Cleaning Procedure  
Column Storage

### **PART VII - REZEX POLYMER-BASED COLUMNS 11-15**

Running Parameters  
Mobile Phase Considerations  
Cleaning Procedure  
Table of Specifications and Operating Parameters

### **PART VIII - POLYSEP-GFC-P COLUMNS 16**

Running Parameters  
Mobile Phase Considerations  
Cleaning Procedure  
Column Storage

### **PART IX - PHENOGE GPC COLUMNS 17, 20-21**

Specifications  
Sample Considerations  
Column Storage

Solvent Switching Considerations  
Solvent Compatibility Chart

### **SOLVENT MISCIBILITY TABLE 18-19**

### **PART X - POLYMERX RP COLUMNS 22**

Specifications  
Running Parameters  
Mobile Phase Considerations

Cleaning Procedure  
Column Storage

### **PART XI - HPLC COLUMN PROTECTION AND PERFORMANCE TESTING 22-29**

Phenex Syringe Filters 22-24  
Membrane Filters Order List Guide 22-23  
Phenex Disposable Centrifugal Filter Units 24  
SecurityGuard Guard Cartridge System 25-29  
HPLC System Test Kit 30  
Column Performance Check Standards 30-33

### **PART XII - SOLID PHASE EXTRACTION (SPE) 34**

### **PART XIII - HPLC ACCESSORIES 35**

### **PHENOMENEX WARRANTY 36**

### **COLUMN PERFORMANCE RECORD 37**

# INTRODUCTION

Every Phenomenex HPLC column is a precision product which, though delicate, will provide excellent performance, reproducibility and column lifetime if cared for properly. The information and recommendations contained in this manual are designed to guide you in the care and use of your column, but should not be considered absolute. Please follow the instructions herein to maximize column performance and lifetime. Should you have any questions, please contact your Phenomenex Technical Representative or local distributor.


## UPON RECEIPT OF THE COLUMN

- **Verify the column you received is the column you ordered**
- **Check the column for physical damage which may have occurred during shipping**
- **Test the column immediately to verify performance and quality**
- **All columns are shipped in the testing solvent, unless otherwise specified**

Each Phenomenex manufactured HPLC column is individually packed and tested to ensure high column quality. Every column is supplied with its Test Chromatogram and a Specification Sheet which indicates column serial number and identity, testing conditions and operating parameters.

The warranty period begins upon receipt of the column. Testing is especially important if the column is to be placed in storage. Test the column using the same conditions in the test chromatogram. Use the formulae in Figure 1 to determine column efficiency and peak asymmetry.

Chromatographic performance depends on the entire system, not just the column. Columns are QC tested using optimum conditions to minimize bandspreading from "extra column effects." Most variations from the Phenomenex test data are due to extra-column effects created by the design of your system (i.e., injector, flow cell, connecting tubing, etc.). If you have any questions regarding your test results or the column quality, or if there are signs of damage, CONTACT PHENOMENEX OR YOUR LOCAL DISTRIBUTOR IMMEDIATELY.



*Formulae for calculating efficiency and peak asymmetry*




## SELECTING THE RIGHT TUBING AND FITTINGS

The tubing and fittings on an HPLC system contribute to system dead volume. If not minimized, dead volume can lead to band broadening and peak degradation. Please use the following guideline to keep system dead volume to a minimum and to help ensure optimum column performance.



### TUBING

The choice of tubing material is based on its chemical resistivity, application and HPLC system considerations (i.e. flow rate, backpressure, etc). Please refer to Tables 1-3 for specifics.

### TUBING DIAMETER

TABLE 1	
<b>High Pressure Tubing:</b>  1/16 in. OD x 0.010 in. ID	
<b>Inlet/Outlet Low Pressure Tubing:</b>	
 1/16 in. OD x 0.030 in. ID	 1/8 in. OD x 0.062 in. ID

### TUBING COMPATIBILITY

TABLE 2		
<b>Stainless Steel</b> (Type 316)		AVOID high concentrations of acids or halogenated salts
<b>PEEK</b> (biocompatible)		AVOID 100 % THF, chlorinated solvents, high concentrations of acids
<b>Titanium</b> (biocompatible)		Compatible with nearly all chemicals


### TUBING APPLICATIONS

TABLE 3		
Tubing ID (Inch)	Column IDs (mm)	Typical Flow Rates (mL/min)
0.002	0.30 (Fused Silica)	0.001 - 0.02
0.005	1.0 (Stainless Steel)	0.02 - 0.1
0.007	2.0 - 4.6	0.2 - 2.0
0.010	3.2 - 7.8	0.5 - 5.0
0.020	10.0 - 21.2	2.0 - 50.0
0.040	21.2 - 100.0	10.0 - 200.0

## FITTINGS

All Phenomenex column end-fittings are female inverted (internal type) with 10-32 type threading:


- The end-fitting can fit any  $\frac{1}{16}$  in. OD tubing (see page 2 for tubing considerations)
- A 10-32 threaded male nut and ferrule or a polymeric fingertight male nut is used to swage or tighten the tubing onto the fitting (see Figure 2)



## INSTALLATION CONSIDERATIONS:

- The shape of the swaged ferrule can differ between manufacturers. For Phenomenex columns, you may use Phenomenex or Valco type ferrules.
- **VERY IMPORTANT:** The seating depth of the nib (Figure 2) for Phenomenex columns is 0.080 in. Tubing **MUST** be seated all the way down into the column end-fitting. Failure to do so will result in having a small mixing chamber at the top or bottom of the column. This will lead to degraded chromatography.

\*Polymeric fingertight fittings are easy to use. They come in one piece, require **NO** tools for attachment and easily conform to the shape of the column end-fitting.



---

## COLUMN INSTALLATION

***IT IS HIGHLY RECOMMENDED THAT YOU READ THIS GUIDE FOR SPECIFIC COLUMN CONSIDERATIONS BEFORE PROCEEDING WITH THE INSTALLATION (PARTS I-X)***

- Flush HPLC pump and line thoroughly with filtered and degassed mobile phase (without any buffers). Make sure there are no air bubbles in the system.
- Connect the column to the injector corresponding to the direction of the flow label (located on the column). Leave the outlet of the column unattached.
- Set pump to flow at 0.1 mL/min (or lowest setting) and increase to normal flow rate over 5 minutes.
- Stop flow when there is a free flow of solvent from the column outlet, wipe the end and attach to the detector
- Equilibrate the column by passing approximately 10-30 column volumes of mobile phase at normal flow rate.
- **For those columns that can be used under reversed-phase or normal phase conditions (i.e., -CN or -NH<sub>2</sub>), flush with 20-30 column volumes of IPA or THF as the intermediate solvent when switching from reversed-phase to normal phase modes, or vice versa.**

## PART I - SILICA-BASED & TWIN™ TECHNOLOGY COLUMNS

### RUNNING PARAMETERS

- Keep backpressures below 3500 psi (245 bar) [maximum 5000 psi] unless otherwise specified in Parts I-X
- Avoid any sudden pressure changes
- If high backpressure is observed, reverse flush the column (do not try this on other manufacturers' columns)
- Use a backpressure regulator if you are experiencing out-gassing problems in the detector cell.
- Maximum operating temperature is 60 °C for all Phenomenex silica-based reversed phase columns.

---

### MOBILE PHASE CONSIDERATIONS

- Use only HPLC grade solvents and water
- Use only highest purity chemicals and reagents
- Filter and degas all mobile phases prior to use
- **Make sure solvents are miscible**

Trace impurities can dramatically degrade HPLC columns. When changing to a different mobile phase, make sure the solvents and/or buffers are miscible (see Table 11). Using solvents that are immiscible with the solvent in the column can permanently damage the column. **Salt and buffer precipitation from the mobile phase can permanently damage the column.** Always check sample solubility and if possible use the mobile phase as the diluent (sample solvent).

---

### STATIONARY PHASE CONSIDERATIONS

- Maintain pH between 2.0 and 8.0\*
- Use presaturator columns and guard columns
- Avoid aldehydes and ketones with amino columns

Silica-based columns are pH sensitive. Low pH ( $\leq 2.0$ ) will hydrolyze the bonded phase (strip off the functional groups) and high pH ( $\geq 8.0$ ) will dissolve the silica. If the mobile phase pH is near 2.0 or 8.0, use a presaturator column.

***\*Consult Phenomenex for columns that have extended pH ranges.***

## BACKPRESSURE AND FLOW RATES

To maximize column life, flow rates should be adjusted to keep pressures below 3500 psi.

TABLE 4

Particle Size $\mu\text{m}$	Internal Diameter(mm)	Typical Flow Rate(mL/min)	Typical Pressure(psi)	
			150 mm*	250 mm*
3	4.6	0.5	985	1640
5	1.0	0.1	1500	2500
5	2.0	0.2	750	1250
5	3.0	0.5	732	1226
5	4.6	1.0	710	1180
5	10.0	5.0	750	1250
10	4.6	2.0	355	590
10	21.2	20.0	170	280

\* column length

Columns can be operated at any flow rate that is consistent with the backpressure limitations described above. Flow rates should be optimized to provide the best efficiency for your sample.

## SCALING UP/SCALING DOWN

Adjusting flow rates for different column internal diameters is straightforward. To keep the retention times constant, the flow rates and loading capacity must be adjusted according to the column's internal diameter. Assuming column length does not change:

$$X = \text{Scale Factor} = \frac{(\text{radius column B})^2}{(\text{radius column A})^2}$$

From a 4.6 mm ID column some approximate scaling factors are:

TABLE 5

Internal Diameter	Scaling Factor
1.0 mm	0.05x
2.0 mm	0.2x
3.0 mm	0.5x
10.0 mm	5x
21.2 mm	21x



HPLC columns running water-free, flammable organic solvents (e.g., normal phase, chiral, GPC) can generate static electricity and should be properly grounded to avoid a potentially dangerous electrical discharge.

## COLUMN STORAGE

- Column storage conditions affect column lifetime
- Never store columns with buffers
- Flush with 5 column volumes of mobile phase without buffer to remove any buffers or salts

Storage Conditions for Silica-Based HPLC Columns:

TABLE 6

Column Type	Storage Solvent
Reversed Phase C18, C12, C8, C4, C2, C1, Phenyl, PFP	65 % Acetonitrile/ 35 % Water
Normal Phase Silica, CN, NH <sub>2</sub> , PAC, Diol Alumina	Isopropanol or Hexane
Ion-Exchange SAX, SCX, WAX, WCX	Methanol*
Size-Exclusion Diol	0.05 % NaN <sub>3</sub> in water or 10 % methanol
HILIC Luna HILIC	80 % Acetonitrile/ 20 % Water

\*Flush column with 50 mL HPLC grade water prior to storage solvent

---

## COLUMN CLEANING PROCEDURES

The following conditions apply to Phenomenex silica-based columns with the exception of chiral columns (see Parts IV and V):

- Before starting any kind of cleaning procedure, make sure your in-column solvent or mobile phase is miscible with the recommended cleaning solvent(s).
- Flow rates should be  $1/5 - 1/2$  of the typical flow rate.

To estimate the column volume, use the following

equation:  $V = \pi r^2 L$

$V$  = column volume in mL

$r$  = column radius in cm

$L$  = column length in cm

### UNBONDED SILICA

Rinse with 10 Column Volumes each of:

- Hexane
- Methylene Chloride
- Isopropanol
- Methylene Chloride
- Mobile Phase

### Water Removal Procedure:

Flush column with 30 mL 2.5 % 2,2-dimethoxy-propane and 2.5 % glacial acetic acid in Hexane

### BONDED NORMAL PHASE (CN, NH<sub>2</sub>, DIOL, PAC)

Rinse with 10 Column Volumes each of:

- Chloroform
- Isopropanol
- Methylene Chloride
- Mobile Phase

*Exception: Luna Amino in reversed phase mode.*

### HILIC

Rinse with 10 Column Volumes each of:

- 95 % Water/5 % Acetonitrile (for buffer removal)
- 95 % 100 mM Ammonium Acetate, pH 5.8/5 % Acetonitrile
- 95 % Water/5 % Acetonitrile
- Mobile Phase

### REVERSED PHASE

**(C18, C12, C8, C4, C5, C2, C1, PHENYL, PFP, CN, NH<sub>2</sub>)**

Rinse with 10 Column Volumes each of:

- 95 % Water/5 % Acetonitrile (for buffer removal)
- THF
- 95 % Acetonitrile/5 % Water
- Mobile Phase

### REVERSED PHASE PROTEIN/PEPTIDE

**(C18, C12, C8, C5, C4, Phenyl)**

Rinse with 20 column volumes of mobile phase with buffer removed. Run gradient (2x):

A) 0.1 % TFA in water

B) 0.1 % TFA in Acetonitrile/Isopropanol (1:2)

25 % B to 100 % B for 30 minutes

Equilibrate with 10 column volumes of mobile phase

Do not store column in TFA

### ION EXCHANGE (SAX, SCX, NH<sub>2</sub>, WAX, WCX)

Rinse with 10 Column Volumes each of:

- 500 mM Phosphate Buffer pH 7
- 10 % Acetic Acid (Aq)
- 5 Column Volumes of Water
- 10 Column Volumes of Phosphate Buffer pH 7
- 5 Column Volumes of Water
- 10 Column Volumes of Methanol
- 10 Column Volumes of Water

### For Protein Removal

Follow the above procedure with this exception:

Substitute 10 Column Volumes of Methanol with 10 Column Volumes of 5 M Urea **or** 5 M Guanidine Thiocyanate



## **GFC/SEC (BioSep SEC\*)** *\*See Part VI for more details*

Rinse with 5 column volumes of 0.1 M Phosphate Buffer pH 3.0. For strongly retained proteins, run the following gradient: 100 % Water to 100 % Acetonitrile to 100 % Water over 60 minutes OR wash with 5 column volumes of SDS or 6 M Guanidine Thiocyanate or 10 % DMSO

---

## **SAMPLE CONSIDERATIONS**

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical to preparative column.

## **PART II - SFC (SUPERCRITICAL FLUID CHROMATOGRAPHY)**

Phenomenex analytical and Axia 'SFC Approved' columns have been leak tested under SFC conditions at pressure far exceeding what may be expected with normal SFC operation.

---

## **RUNNING PARAMETERS**



- Backpressure Limitation: 3500 psi
- Flow Rate Limitations: Flow rate is to be controlled so that pressure limit of 3500 psi is not exceeded
- pH Limitations: Dictated by the media packed in the column

---

## **EQUILIBRATING COLUMN**

SFC column stationary phases have a polar surface and may be shipped under reversed phase or normal phase conditions. Flush all columns with 10-30 column volumes of Methanol/CO<sub>2</sub> as intermediate solvent between CO<sub>2</sub> and column shipping conditions. Be aware of backpressure settings.

Equilibrate column to starting conditions with 10 column volumes of mobile phase.

---

## **MOBILE PHASE CONSIDERATIONS**

- Use only HPLC grade solvent modifiers
- Use only highest purity chemicals and reagents
- Filter and degas all mobile phases prior to use

---

## **CLEANING PROCEDURE**

- Under extreme conditions the column can be flushed with 50/50 Acetonitrile/Isopropyl Alcohol followed by 100 % Isopropyl Alcohol. Maintain backpressure below limits.
- Re- Equilibrate column to starting conditions with 10 column volumes of mobile phase

---

## **COLUMN STORAGE**

- Completely remove all buffers, acids, bases or other mobile phase additives to prevent damage to media
- Flush with at least 10 column volumes of Methanol after the last sample is purified
- Store column with end plugs firmly seated in endfittings to ensure storage solvent does not evaporate

## PART III - AXIA PACKED PREPARATIVE COLUMN

### RUNNING PARAMETERS

- Backpressure Limitation: 3500 psi
- Flow Rate Limitations: Determined by the viscosity of mobile phase; flow rates to be controlled so that backpressure limit of 3500 psi is not exceeded
- pH Limitations: Dictated by the media packed in the column

---

### MOBILE PHASE CONSIDERATIONS

- Use only HPLC grade solvents and water
- Use only highest purity chemicals and reagents
- Filter and degas all mobile phases prior to use

---

### CLEANING PROCEDURE

(for Axia reversed phase columns)

- Reverse flush the column with 10 column volumes of HPLC grade water and then 10 column volumes of 100 % organic solvent
- Under extreme conditions, the column can be flushed with 10 column volumes of 100 % THF (or IPA) followed by 100 % methylene chloride
- After cleaning, wash with 100 % THF (or IPA) and 50:50 Acetonitrile/Water, prior to equilibrating with the starting mobile phase

---

### COLUMN STORAGE

- Completely remove all buffers, acids, bases, or other mobile phase additives to prevent physical damage to the media
- Flush with at least 10 column volumes of 50:50 Acetonitrile/Water after the last sample is purified
- Store with column end plugs placed back in the end-fittings to ensure that the packing media does not dry out

*For additional information, consult the Care and Use of Axia Packed Preparative HPLC Columns, included with each Axia column purchased.*

## PART IV - LUX CHIRAL COLUMNS

### RUNNING PARAMETERS

#### OPERATING BACKPRESSURE

The mobile phase flow rate should be set such that the column backpressure stays below 300 bar (4300 psi). This maximum backpressure should not be exceeded for long periods of time.

#### OPERATING TEMPERATURES

With standard mobile phases (such as alkane/alcohol) the column can be used in the temperature range 0-50 °C.

---

### MOBILE PHASE CONSIDERATIONS

#### MOBILE PHASE COMPATIBILITY

Lux columns can be used with normal phase (alkane/alcohol), reversed phase (aqueous methanol, aqueous acetonitrile or appropriate buffer/methanol or buffer/acetonitrile mixtures), as well as with pure polar organic solvents (low molecular weight alcohols, acetonitrile or their mixtures).

#### SOLVENT SWITCHING

An appropriate column washing procedure must be applied when changing from one mobile phase to another. The miscibility of the different mobile phase components must be carefully considered for this wash. To safely transfer a column from hexane

to methanol (or acetonitrile) or from methanol (or acetonitrile) to hexane, use 100 % 2-propanol as transition solvent at a flow rate of 0.2-0.5 mL/min. Ten column volumes of 2-propanol (i.e. 25 mL for a 250 x 4.6 mm ID column or 15 mL for a 150 x 4.6 mm ID column) are sufficient for completely removing the old mobile phase. In addition, when the buffer salt additive of the RP mobile phase is insoluble in 2-propanol, flush the column briefly with water before switching to a buffered mobile phase.

Switching between elution modes – NP to PO or RP and back to NP – is possible, but changes in resolution and retention times may be observed with some compounds. We recommend the use of dedicated Lux columns to reversed phase operation hence avoiding the need of converting columns used in normal phase elution mode to reversed phase or vice versa.

## **USE OF MOBILE PHASE MODIFIERS**

For basic samples or acidic chiral compounds, it may be necessary to use an appropriate mobile phase modifier in order to achieve chiral resolution and to insure proper peak shapes. Diethylamine, ethanolamine and butyl amine in the concentration range 0.1-0.5 % can be used with basic analytes, while trifluoroacetic or acetic acid (0.1-0.5 %; typically 0.1-0.2 %) with acidic analytes. Mixtures of basic and acidic mobile phase additives are acceptable (e.g. diethylamine acetate or trifluoroacetate). Lux columns will deliver consistent results when operated with mobile phases containing additives at the concentration levels specified above. However, limited decrease in column efficiency may occur when a column is used in combination with these additives. Therefore, we advise to dedicate columns to mobile phases containing basic additives.

## **MOBILE PHASE RESTRICTIONS**

Lux chiral stationary phases are prepared by coating silica with various polysaccharide derivatives. Therefore, any solvent dissolving the polysaccharide derivative (such as tetrahydrofuran, acetone, chlorinated hydrocarbons, ethylacetate, dimethylsulfoxide, dimethylformamide, N-methylformamide, etc.) must be avoided even in trace amounts (e.g. even as sample solvent).

---

## **EXTENDING LIFETIME AND RECONDITIONING**

Phenomenex recommends the use of SecurityGuard™ guard cartridges to extend the lifetime of your column, especially with samples extracted from complex matrixes. Ideally, samples must be completely dissolved in the mobile phase or filtered through a syringe filter of approximately 0.45 µm porosity.

---

## **COLUMN STORAGE**

- Column storage for a longer period of time is recommended in n-hexane/2-propanol (9:1, v/v).
- Columns used in reversed phase conditions should be first flushed with water (whenever a buffer salt was used as RP mobile phase additive) and then with methanol (or with methanol only when no salt was used). The column can be stored in methanol.

## PART V - CHIREX CHIRAL COLUMNS

### RUNNING PARAMETERS

- Temperature must not exceed 50 °C
- Column pressure must not exceed 3000 psi
- Maintain flow rate between 0.5-2.0 mL/min for 4.6 mm ID columns

---

### MOBILE PHASE CONSIDERATIONS

- Dedicate column to reversed or normal phase solvents
- pH range: 2.5 to 7.5
- Use only HPLC grade solvents
- Use only highest purity chemicals and reagents
- Filter and degas all mobile phases prior to use
- Make sure solvents are miscible (see pp. 18-19)

*Most CHIREX Chiral columns use a Type I or brush type Chiral stationary phase (CSP I). Normal phase systems usually provide better selectivity than reversed phase systems. SEE COLUMN INSERT FOR FURTHER INFORMATION ABOUT SPECIFIC CHIREX COLUMNS (included with each column)*

---

### SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical to preparative column.

## PART VI - BIOSEP-SEC-S COLUMNS

### RUNNING PARAMETERS

- Maximum flow rate is a function of pressure
- Column pressure must not exceed 1500 psi
- Maximum temperature: 50 °C for 316 stainless steel

---

### MOBILE PHASE CONSIDERATIONS

- pH range: 2.5 - 7.5
- Maximum organic modifier: Up to 100 % CH<sub>3</sub>CN. Start with 100 % H<sub>2</sub>O, linear gradient to 100 % CH<sub>3</sub>CN over 50 min. Up to 90 % CH<sub>3</sub>CN, 10 % DMSO or 500 mM β-mercaptoethanol.
- Maximum salt concentration: 1 M
- Filter and degas all mobile phases prior to use.

---

### SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical to preparative column.

---

### CLEANING PROCEDURE

- General protein removal: wash with 30 mL of 0.1 M NaH<sub>2</sub>PO<sub>4</sub>, pH 3.0
- Hydrophobic protein removal: use Acetonitrile gradient
- Strongly adsorbed proteins: wash with 30 mL of 0.5 % SDS or 6 M guanidine thiocyanate, or 10 % DMSO

---

### COLUMN STORAGE

- Overnight storage: run mobile phase at 0.2 mL/min
- Prolonged storage: use 0.05 % sodium azide in water or 10 % methanol in water

# PART VII - REZEX POLYMER-BASED COLUMN

## RUNNING PARAMETERS

- Columns must be run at elevated temperatures
- (60-85 °C) except Rezex ROA and RHM for most applications
- Column pressure for 8 % cross-linked material must not exceed 1,000 psi; must not exceed 300 psi for 4 % cross-linked material
- Clean and reverse flush column regularly with HPLC grade water

**Important:** Never exceed maximum pressure limitations. This will cause irreversible damage to the column.

---

## MOBILE PHASE CONSIDERATIONS

- Filter and degas all mobile phases prior to use
- Do not exceed 10 % Organic, IPA, EtOH
- Store columns in HPLC grade water

Rezex utilizes a sulfonated polystyrene resin which is very rugged and resistant to chemical attack. However, the material is pressure sensitive and must be cared for properly.

### START UP

Turn on column heating unit to 60 - 85 °C and start the mobile phase at 0.1 mL/min. Make sure the pressure remains below 400 psi for 8 % cross-linked material; below 200 psi for 4 % cross-linked material. As the temperature reaches working condition, increase flow rate to the specified level. (See Rezex Operating Parameters)

### SHUT DOWN

Overnight: Lower flow rate to 0.1 mL/min. Leave system on and continue heating.

Long Term: Store columns in 100 % water. Turn off pump and allow the system to cool. Replace the end plugs and tightly cap the column.

---

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical column.

---

## CLEANING PROCEDURE

Before utilizing any cleaning procedure outlined in the Tables on pages 13 and 15, first try to clean your Rezex column as follows:

Remove the guard column and reverse the direction of flow on the analytical column. Run 100 % HPLC grade water through the column as follows:

TABLE 7

Rezex Column	Flow (mL/min)	Temp. (°C)
RPM, RCM, RHM	0.6	85
RCU	0.2	85
RSO and RNO	0.2	75
RNM and RAM	0.4	75
ROA	0.6	85

Run the column under these conditions for a minimum of 12 hours. After completing the cleaning procedure, return the column to the original direction of flow and equilibrate for analysis.

If this procedure is not effective in cleaning the column, proceed to the specified procedures outlined in Tables 8 and 9.

# PART VII - REZEX POLYMER-BASED COLUMNS (cont'd)


## SPECIFICATIONS AND OPERATING PARAMETERS

<b>Table 8</b>	<b>RCM Monosaccharide</b>	<b>RSO Oligosaccharide</b>	<b>RNO Oligosaccharide</b>	<b>RNM Carbohydrate</b>	<b>RAM Carbohydrate</b>
<b>Part Number</b>	00H-0130-K0	00P-0133-N0	00P-0137-N0	00H-0136-K0	00H-0131-K0
<b>Ionic Form</b>	Calcium	Silver	Sodium	Sodium	Silver
<b>Standard Dimensions</b>	300 x 7.8 mm	200 x 10 mm	200 x 10 mm	300 x 7.8 mm	300 x 7.8 mm
<b>Matrix</b>		Sulfonated Styrene Divinyl Benzene			
<b>Cross Linking</b>	8 %	4 %	4 %	8 %	8 %
<b>Particle Size (µm)</b>	8	12	12	8	8
<b>Min. Efficiency (p/m) based on last peak</b>	35,000	N/A	N/A	30,000	35,000
<b>Typical Pressure (psi @ Max Flow Rate)</b>	400	200	200	400	400
<b>Max. Pressure (psi @ Max Flow Rate)</b>	1,000	300	300	1,000	1,000
<b>Max. Flow Rate (mL/min)*</b>	1.0	0.3	0.3	1.0	1.0
<b>Max. Temperature (°C)</b>	85	85	85	85	85
<b>Typical Mobile Phase</b>	Water	Water	Water	Water	Water
<b>pH Range</b>	Neutral	Neutral	Neutral	Neutral	Neutral
<b>Guard Column Part No.</b>	03B-0130-K0	03R-0133-N0	03R-0137-N0	03B-0136-K0	03B-0131-K0

\* Make sure the maximum pressure is not exceeded

# COLUMNS FOR CARBOHYDRATE AND ORGANIC ACID ANALYSIS

**Table 8 (continued)**

<b>Cleaning, Regeneration and Storage</b>		<b>RCM Monosaccharide</b>	<b>RSO Oligosaccharide</b>	<b>RNO Oligosaccharide</b>	<b>RNM Carbohydrate</b>	<b>RAM Carbohydrate</b>
Organic Modifiers (Max)		10 % Methanol, IPA, EtOH, Acetonitrile				
Inorganic Modifiers (Max)		5 % CaSO <sub>4</sub> , Ca(NO <sub>3</sub> ) <sub>2</sub> , CaCl <sub>2</sub>	5 % Silver Nitrate	5 % Sodium Salts	5 % Sodium Salts	2 % Silver Nitrate
Avoid		Acids, Bases, Non-Calcium Salts or Metal Ions, >10 % Organic	Acids, Bases, Non-Silver Salts/Metal Ions, >5 % Organic	Acids, Bases, Non-Sodium Salts/Metal Ions, >5 % Organic	Acids, Bases, Non-Sodium Salts/Metal Ions, >10 % Organic	Acids, Bases, Non-Silver Salts/Metal Ions, >10 % Organic
Cleaning Solvent		100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate (mL/min)		0.6	0.2	0.2	0.4	0.4
Temperature (°C)		85	75	75	75	75
Duration (hrs)		12	12	12	12	12
Regeneration Solvent		0.1 M Ca(NO <sub>3</sub> ) <sub>2</sub>	0.1 M AgNO <sub>3</sub>	0.1 M NaNO <sub>3</sub>	0.1 M NaNO <sub>3</sub>	0.1 M AgNO <sub>3</sub>
Flow Rate (mL/min)		0.2	0.1	0.2	0.2	0.2
Temperature (°C)		85	85	85	85	85
Duration (hrs)		4-16	4-16	4-16	4-16	4-16
Ship/Storage Solvent		Water	Water	Water	Water	Water

# PART VII - REZEX POLYMER-BASED COLUMNS (cont'd)

## SPECIFICATIONS AND OPERATING PARAMETERS

<b>Table 9</b>	<b>RPM Monosaccharide</b>	<b>RHM Monosaccharide</b>	<b>ROA Organic Acid</b>	<b>RFQ Fast Acid</b>	<b>RCU Sugar Alcohols</b>
<b>Part Number</b>	00H-0135-K0	00H-0132-K0	00H-0138-K0	00D-0223-K0	00G-0130-D0
Ionic Form	Lead	Hydrogen	Hydrogen	Hydrogen	Calcium
Standard Dimensions	300 x 7.8 mm	300 x 7.8 mm	300 x 7.8 mm	100 x 7.8 mm	250 x 4.0 mm
Matrix	Sulfonated Styrene Divinyl Benzene				
Cross Linking	8 %	8 %	8 %	8 %	8 %
Particle Size (µm)	8	8	8	8	8
Min. Efficiency (p/m) (based on last peak)	35,000	35,000	50,000 (Acetic Acid)	30,000	12,000
Typical Pressure (psi @ Max Flow Rate)	400	400	400	400	400
Max. Pressure (psi @ Max Flow Rate)	1,000	1,000	1,000	1,000	1,000
Max. Flow Rate (mL/min)*	1.0	1.0	1.0	1.0	0.5
Max. Temperature (°C)	85	85	85	85	85
Typical Mobile Phase	Water	Water	0.005N H <sub>2</sub> SO <sub>4</sub>	0.005N H <sub>2</sub> SO <sub>4</sub>	Water
pH Range	Neutral	1-8	1-8	1-8	Neutral
Guard Column Part No.	03B-0135-K0	03B-0132-K0	03B-0138-K0	03B-0223-K0	03A-0130-D0

\* Make sure the maximum pressure is not exceeded



# COLUMNS FOR CARBOHYDRATE AND ORGANIC ACID ANALYSIS

**Table 9 (continued)**

<b>Cleaning, Regeneration and Storage</b>		<b>RPM Monosaccharide</b>	<b>RHM Monosaccharide</b>	<b>ROA Organic Acid</b>	<b>RFQ Fast Acid</b>	<b>RCU Sugar Alcohols</b>
Organic Modifiers (Max)		10 % Methanol, IPA, EtOH, Acetonitrile				
Inorganic Modifiers (Max)		5 % Lead Nitrate	5 % HNO <sub>3</sub> , H <sub>3</sub> PO <sub>4</sub>	5 % HNO <sub>3</sub> , H <sub>3</sub> PO <sub>4</sub>	5 % HNO <sub>3</sub> , H <sub>3</sub> PO <sub>4</sub>	5 % CaSO <sub>4</sub> , Ca(NO <sub>3</sub> ) <sub>2</sub> , CaCl <sub>2</sub>
Avoid		Acids, Bases, Non-Lead Salts/Metal Ions, >10 % Organic	Acids, Bases, Salts, Metal Ions, >10 % Organic	Acids, Bases, Salts, Metal Ions, pH > 3, >10 % Organic	Acids, Bases, Salts, Metal Ions, pH > 3, >10 % Organic	Acids, Bases, Non-Calcium Salts or Metal Ions, >10 % Organic
Cleaning Solvent		100 % Water	100 % Water	100 % Water	100 % Water	100 % Water
Flow Rate (mL/min)		0.6	0.6	0.6	0.6	0.2
Temperature (°C)		85	85	85	85	85
Duration (hrs)		12	12	12	12	12
Regeneration Solvent		0.1 M Pb(NO <sub>3</sub> ) <sub>2</sub>	0.025 M H <sub>2</sub> SO <sub>4</sub>	0.025 M H <sub>2</sub> SO <sub>4</sub>	0.025 M H <sub>2</sub> SO <sub>4</sub>	0.1 M Ca (NO <sub>3</sub> ) <sub>2</sub>
Flow Rate (mL/min)		0.2	0.2	0.2	0.2	0.2
Temperature (°C)		85	85	85	85	85
Duration (hrs)		4-16	4-16	4-16	4-16	4-16
Ship/Storage Solvent		Water	Water	0.005 N H <sub>2</sub> SO <sub>4</sub>	0.005 N H <sub>2</sub> SO <sub>4</sub>	Water

See p. 11 for general care and usage of Rezex columns.

# PART VIII - POLYSEP-GFC-P COLUMNS

## RUNNING PARAMETERS

- Column pressure must not exceed 650 psi
- Do not exceed 60 °C

---

## MOBILE PHASE CONSIDERATIONS

- pH range: 3 - 12
- Maximum salt concentration: 0.5 M
- Organic Modifier capacity:

	<i>POLYSEP PHASE</i>						
	<b>1000</b>	<b>2000</b>	<b>3000</b>	<b>4000</b>	<b>5000</b>	<b>6000</b>	<b>Linear</b>
<b>Methanol</b>	20 %	95 %	70 %	70 %	70 %	70 %	70 %
<b>Acetonitrile</b>	20 %	70 %	70 %	70 %	70 %	70 %	70 %

---

## CLEANING PROCEDURE

0.5 % SDS or 6 M guanidine thiocyanate. All PolySep columns except for PolySep 1000 may also be cleaned with 50 % acetonitrile. Make sure not to exceed a maximum pressure of 650 psi when cleaning.

---

## COLUMN STORAGE

- Overnight storage: run buffer at low flow rate (0.2 mL/min or less)
- Prolonged storage: store in 0.05 % sodium azide in water or 10 % methanol in water

---

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the column. Use of a guard column is highly recommended to prolong the life of your analytical to preparative column.

# PART IX - PHENOGEL GPC COLUMNS

## SPECIFICATIONS

<b>Matrix:</b>	Styrene-Divinyl Benzene Copolymer
<b>Particle Size:</b>	5, 10, 20 $\mu\text{m}$
<b>Porosities:</b>	50 $\text{\AA}$ to $10^6 \text{\AA}$ , and mixed beds
<b>Typical Pressure:</b>	<b>5 <math>\mu\text{m}</math>:</b> 300 psi <b>10 <math>\mu\text{m}</math>:</b> 200 psi
<b>Maximum Pressure:</b>	650 psi
<b>Maximum Temperature:</b>	140 $^{\circ}\text{C}$ (205 $^{\circ}\text{C}$ for UT)
<b>Minimum Efficiency*:</b>	<b>5 <math>\mu\text{m}</math>:</b> 45,000 P/m** <b>10 <math>\mu\text{m}</math>:</b> 35,000 P/m**
<b>Typical Flow Rates:</b>	<b>4.6 mm ID:</b> 0.35 mL/min <b>7.8 mm ID:</b> 1.0 mL/min <b>21.2 mm ID:</b> 7.0 mL/min
<b>End Fittings:</b>	Valco Compatible

\*Tested in THF

\*\* For 300 x 7.8 mm ID columns

## SAMPLE CONSIDERATIONS

Always prefilter samples with Phenex™ syringe filters to avoid particulate contaminants which may clog the GPC column. Use of a GUARD COLUMN is highly recommended to prolong the life of your analytical or preparative column. For optimal results, use the chart below to determine sample concentrations and injection volumes.

TABLE 10

Molecular Weight	Concentration (w/v)	Max Injection Volume
< 50 K	0.5 %	100 $\mu\text{L}$
50-600 K	0.25 %	100 $\mu\text{L}$
600-3000 K	0.05 %	100 $\mu\text{L}$
>3000 K	0.01 %	20 $\mu\text{L}$

Continued on p. 20

# SOLVENT MISCIBILITY TABLE

TABLE 11

## Solvent Polarity Chart

Relative Polarity	Compound Formula	Group	Representative Solvent Compounds
<b>NONPOLAR</b>	R - H	Alkanes	Petroleum ethers, ligroin, hexanes
	Ar - H	Aromatics	Toluene, benzene
	R - O - R	Ethers	Diethyl ether
	R - X	Alkyl halides	Tetrachloromethane, chloroform
	R - COOR	Esters	Ethyl acetate
	R - CO - R	Aldehydes and ketones	Acetone, methyl ethyl ketone (MEK)
	R - NH <sub>2</sub>	Amines	Pyridine, triethylamine
	R - OH	Alcohols	Methanol, ethanol, isopropanol, butanol
	R - COHN <sub>2</sub>	Amides	Dimethylformamide
	R - COOH	Carboxylic acids	Ethanoic acid
<b>POLAR</b>	H - OH	Water	Water

Increasing Polarity

Solvent	Polarity Index	Refractive Index @ 20 °C	UV (nm) Cutoff @ 1 AU	Boiling Point (°C)	Viscosity (cPoise)	Solubility in Water (% w/w)
Acetic Acid	6.2	1.372	230	118	1.26	100
Acetone	5.1	1.359	330	56	0.32	100
Acetonitrile	5.8	1.344	190	82	0.37	100
Benzene	2.7	1.501	280	80	0.65	0.18
Butyl Acetate	4.0	1.394	254	125	0.73	0.43
n-Butanol	3.9	1.399	215	118	2.98	7.81
Carbon tetrachloride	1.6	1.466	263	77	0.97	0.08
Chloroform	4.1	1.446	245	61	0.57	0.815
Cyclohexane	0.2	1.426	200	81	1.00	0.01
1,2-Dichloroethane <sup>1</sup>	3.5	1.444	225	84	0.79	0.81
Dichloromethane <sup>2</sup>	3.1	1.424	235	41	0.44	1.6
Dimethylformamide	6.4	1.431	268	155	0.92	100
Dimethyl sulfoxide <sup>3</sup>	7.2	1.478	268	189	2.00	100
Dioxane	4.8	1.422	215	101	1.54	100
Ethyl Acetate	4.4	1.372	260	77	0.45	8.7
Ethanol	5.2	1.360	210	78	1.20	100
di-Ethyl Ether	2.8	1.353	220	35	0.32	6.89
Heptane	0.0	1.387	200	98	0.39	0.0003

Xylene		2.5	1500	290	0.61	0.018
Water		9.0	1333	200	1.00	100
Trichloroethylene		1.0	1477	273	0.57	0.11
Toluene		2.4	1496	285	0.59	0.051
Tetrahydrofuran		4.0	1407	215	0.55	100
di-iso-Propyl Ether		2.2	1368	220	0.37	
iso-Propanol <sup>6</sup>		3.9	1377	210	2.30	100
n-Propanol		4.0	1384	210	2.27	100
Pentane		0.0	1358	200	0.23	0.004
Methyl Ethyl Ketone <sup>5</sup>		4.7	1379	329	0.45	24
Methyl-t-Butyl Ether <sup>4</sup>		2.5	1369	210	0.27	4.8
Methanol		5.1	1329	205	0.60	100
Hexane		0.0	1375	200	0.33	0.001

& **Miscible**  
**Immiscible\***

### SYNONYM TABLE

<sup>1</sup> Ethylene Chloride	<sup>4</sup> tert-Butyl Methyl Ether
<sup>2</sup> Methylene Chloride	<sup>5</sup> 2-Butanone
<sup>3</sup> Methyl Sulfoxide	<sup>6</sup> 2-Propanol

\*Immiscible means that in some proportions two phases will be produced.

## COLUMN STORAGE

Solvents such as THF (stabilized THF only), Chloroform, Methylene Chloride, and Toluene are commonly used for column storage. Be sure to follow solvent switching instructions (see below) if using solvents other than THF. Storage solvents that remain liquefied at ambient temperatures and are not oxidizing can be used for storage.


**BE SURE THAT ANY COLUMN THAT IS NOT USED IS CAPPED TIGHTLY WITH END-PLUGS TO AVOID EVAPORATION OF SOLVENTS FROM COLUMN. COLUMN DESICCATION IS THE MOST COMMON SOURCE OF COLUMN FAILURE.**

## SOLVENT SWITCHING CONSIDERATIONS FOR NON-AQUEOUS GPC COLUMNS

Phenogel columns are rugged and exhibit wide solvent compatibility. Different solvents, however, produce different swell characteristics (Table 12). Improper solvent switches can result in a void. For this reason, we recommend that you dedicate columns to specific solvents.

If you need to switch solvents, it is VERY IMPORTANT to take the following into consideration:

1. Reduce flow rate to 0.2 mL/min.
2. Backpressure must NEVER exceed 650 psi.
3. Always check solvent miscibility in a beaker or follow the solvent miscibility table on page 18-19 before proceeding with ANY solvent switch.
4. Compare the swell characteristics of solvent 1 (old solvent) to solvent 2 (new solvent) and use the following guideline:
  - If the solvent 1 and solvent 2 belong to the same swell category (Table 12), check the solvent miscibility and proceed with the switch.
  - If solvent 1 and solvent 2 belong to successive swell categories as indicated by the arrows on Table 12, check the miscibility and proceed with the switch.
  - If solvent 1 and solvent 2 DO NOT belong to the same OR successive swell categories, switch to an intermediate solvent FIRST, as indicated by the arrows on Table 12.



# SOLVENT COMPATIBILITY CHART FOR PHENOGEL GPC COLUMNS

Mobile Phase Solvent	PHENOGEL PORE SIZE								Suggested Operating Temperature
	50 Å	100 Å	500 Å	10 <sup>3</sup> Å	10 <sup>4</sup> Å	10 <sup>5</sup> Å	10 <sup>6</sup> Å	Linear & Mixed	
Acetone	Y	Y	Y	Y	Y	Y	Y	Y	
Benzene	Y	Y	Y	Y	Y	Y	Y	Y	
Carbon Tetrachloride	Y	Y	Y	Y	Y	Y	Y	Y	
Chloroform	Y	Y	Y	Y	Y	Y	Y	Y	
30 % HFIP/chloroform	Y	Y	Y	Y	Y	Y	Y	Y	
Diethyl Ether	Y	Y	Y	Y	Y	Y	Y	Y	
Dimethylacetamide (DMAC)	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Dimethylformamide (DMF)	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Dioxane	Y	Y	Y	Y	Y	Y	Y	Y	
DMSO	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Ethyl Acetate	Y	Y	Y	Y	Y	Y	Y	Y	
Hexafluoroisopropanol (HFIP)	Y	Y	Y	Y	Y	Y	Y	Y	

\*Not recommended on 5 µm 50 Å columns.

N = Not Compatible Y = Compatible

Table 13

Mobile Phase Solvent	PHENOGEL PORE SIZE								Suggested Operating Temperature
	50 Å	100 Å	500 Å	10 <sup>3</sup> Å	10 <sup>4</sup> Å	10 <sup>5</sup> Å	10 <sup>6</sup> Å	Linear & Mixed	
Hexane	Y	Y	Y	Y	Y	Y	Y	Y	
m-Cresol	Y*	Y	Y	Y	Y	Y	Y	Y	100 °C
Methyl Ethyl Ketone	Y	Y	Y	Y	Y	Y	Y	Y	
Methylene Chloride	Y	Y	Y	Y	Y	Y	Y	Y	
o-Chlorophenol	Y*	Y	Y	Y	Y	Y	Y	Y	100 °C
o-Dichlorobenzene	Y*	Y	Y	Y	Y	Y	Y	Y	135 °C
Quinolin	Y*	Y	Y	Y	Y	Y	Y	Y	60 °C
Tetrahydrofuran	Y	Y	Y	Y	Y	Y	Y	Y	
Toluene	Y	Y	Y	Y	Y	Y	Y	Y	
Trichlorobenzene	Y*	Y	Y	Y	Y	Y	Y	Y	135 °C
Water	N	N	N	N	N	N	N	N	
Xylene	Y	Y	Y	Y	Y	Y	Y	Y	

## PART X - POLYMERX RP COLUMNS

### SPECIFICATIONS

- Matrix: Polystyrene Divinylbenzene (PSDVB)
- Particle Size: 3, 5, 7, 10  $\mu\text{m}$
- Pore Size: 100  $\text{\AA}$

### RUNNING PARAMETERS

- Maximum temperature: 60  $^{\circ}\text{C}$
- Maximum pressure: 2500 psi

### MOBILE PHASE CONSIDERATIONS

- pH range: 0 - 14
- Avoid buffer strength  $> 0.5 \text{ N}$

### CLEANING PROCEDURE

- 100 % Water to 100 % Acetonitrile. Repeat 3 times.

### COLUMN STORAGE

- 75:25 Acetonitrile / Water

## PART XI - HPLC COLUMN PROTECTION & PERFORMANCE TESTING

- Maximize the life of your valuable HPLC Column
- Reduce system wear and tear
- Save time and money

### PHENEX™ SYRINGE FILTERS

- Increase column lifetime (save money!)
- Ensure more accurate, consistent results
- Eliminate damaging microparticulates

Particulates can damage expensive equipment, valves, columns and pumps. They can also lead to erratic analytical results. Pre-filtering samples prior to analysis is critical in preventing column and frit blockage, undue wear on valve seals, and abnormally high operating pressures.

TABLE 14

Sample or Mobile Phase Volume (mL)	Filter Membrane (diameter, mm)	Format
$\leq 2$	4	Syringe filter
2 to 10	15-17	Syringe filter
10 to 100	25-28	Syringe filter
$> 100$	47	Membrane disk
$> 1000$	90	Membrane disk

### MEMBRANE FILTERS ORDER LIST GUIDE

#### REGENERATED CELLULOSE (RC)

As a universal hydrophilic membrane, RC is widely used in chromatography for the clarification of aqueous samples and solvents. Due to its ultra-low binding capabilities, RC membranes are an excellent choice for proteins, peptides and other biomolecules.

#### POLYTETRAFLUOROETHYLENE (PTFE, TEFLON®)

PTFE is an inherently hydrophobic membrane, excellent for filtration of organic-based, highly acidic or basic samples and solvents. Widely used in chromatography, it is especially well suited for the clarification of non-aqueous samples. Although this membrane is hydrophobic, it can be made hydrophilic by



wetting the membrane with alcohol and then flushing with deionized water.

## POLYETHERSULFONE (PES)

Polyethersulfone, a hydrophilic membrane with fast flow, high-throughput characteristics, with ultra-low protein binding. It is ideally suited for use in life sciences applications. The PES membrane offers better chemical resistance than cellulose acetate. Recommended for filtering critical biological sampling, tissue culture media, additives, and buffers.

## NYLON (NY)

Nylon has inherent hydrophilic characteristics and works well for filtration of many aqueous and mixed-organic samples. Nylon exhibits a high non-specific affinity for proteins. Phenomenex recommends Phenex-RC (Regenerated Cellulose) filters for application requiring low non-specific adsorption of proteins.

## CELLULOSE ACETATE (CA)

Cellulose Acetate membranes exhibit ultra-low protein binding and are broadly used in the filtration of biological samples. In combination with a glass pre-filter (Phenex-GF/CA), this membrane is excellent for filtration of tissue culture media, general biological sample filtration and clarification.

## GLASS FIBER (GF)

Glass Fiber filters are made of inert borosilicate glass and have a nominal 1.2  $\mu\text{m}$  pore size. They are commonly used with highly viscous samples or samples containing high concentrations of particulate matter (e.g., food analysis, biological samples, soil samples, fermentation broth samples, removal of yeasts, molds, etc.). Glass Fiber filters can be used alone or in conjunction with other Phenex filter membranes such as the 0.45  $\mu\text{m}$  pore Phenex-RC filter to reduce clogging of the membrane and optimize flow.

## ORDERING INFORMATION

Part No.	Pore Size ( $\mu\text{m}$ )	Phenex Membrane	Housing
<b>4 mm Diameter (500/pk)</b>			
AF0-3103-52	0.45	RC	PP
AF0-3102-52	0.45	PTFE <sup>6</sup>	PP
AF3-3107-52	0.45	NY	PP
AF0-3203-52	0.20	RC	PP
AF0-3202-52	0.20	PTFE <sup>6</sup>	PP
AF3-3207-52	0.20	NY	PP
<b>15–17 mm Diameter (500/pk)</b>			
AF0-2103-52	0.45	RC	PP
AF0-2102-52	0.45	PTFE <sup>6</sup>	PP
AF0-2107-52	0.45	NY	PP
AF0-2203-52	0.20	RC	PP
AF0-2202-52	0.20	PTFE <sup>6</sup>	PP
AF0-2207-52	0.20	NY	PP
<b>25–28 mm Diameter (500/pk)</b>			
AF0-8103-52 <sup>5</sup>	0.45	RC	PP
AF0-8108-52 <sup>7</sup>	0.45	PES <sup>3</sup>	PP
AF0-1102-52	0.45	PTFE <sup>6</sup>	PP
AF0-1107-52	0.45	NY	PP
AF0-8B09-52 <sup>7</sup>	0.45	GF/CA <sup>2,3,4</sup>	MBS
AF0-8203-52 <sup>5</sup>	0.20	RC	PP
AF0-8208-52 <sup>7</sup>	0.20	PES <sup>3</sup>	PP
AF0-1202-52	0.20	PTFE <sup>6</sup>	PP
AF0-1207-52	0.20	NY	PP
AF0-8A09-52 <sup>7</sup>	0.20	GF/CA <sup>2,3,4,7</sup>	MBS
AF0-8515-52 <sup>7</sup>	1.20	GF <sup>2,3</sup>	MBS


*Housing is made of medical-grade polypropylene (PP), unless otherwise indicated. Above syringe filters are non-sterile.*

- 17 mm diameter.
- Glass fiber filters are 28 mm diameter and made of borosilicate. They will remove 90 % of all particles >1.2  $\mu\text{m}$ .
- Housing material is methacrylate butadiene styrene (MBS) polymerisate. Also known as Cryolite™.
- Cellulose acetate is surfactant-free.
- 26 mm diameter.
- Hydrophobic membrane. Can be made hydrophilic by pre-wetting with IPA.
- 28 mm diameter.

---

## PHENEX™ DISPOSABLE CENTRIFUGAL FILTER UNITS

- Convenient filtration of multiple HPLC and GC samples
- High recovery for small samples
- Nylon, Cellulose Acetate, and PTFE (Teflon®) membrane materials



Centrifugal force drives the sample through the filter quickly without effort on the part of the chemist. No cleaning of syringes is required between samples. The receiver tube serves as a container for the filtered sample and can be retained as long as desired.

---

## ORDERING INFORMATION

Part No.	Pore Size (µm)	Volumes (mL) Sample/Receiver	Membrane Non-Sterile	Unit
AF0-0438	0.2	2.0 / 5.0	Nylon	25/pk
AF0-0439	0.45	2.0 / 5.0	Nylon	25/pk
AF0-0440	0.2	2.0 / 5.0	PTFE	25/pk
AF0-0441	0.45	2.0 / 5.0	PTFE	25/pk
AF0-8353	0.2	2.0 / 5.0	CA	25/pk
AF0-8354	0.45	2.0 / 5.0	CA	25/pk

*Above centrifugal filters are non-sterile.*

# GUARD CARTRIDGE SYSTEM



SecurityGuard provides a great balance of convenience, column protection capability and value. If you've ever used another guard cartridge system or conventional guard column, you will be pleasantly surprised when you see how practical and effective SecurityGuard really is. This highly advanced, patented design offers several unique features up to now not available.

## CONVENIENCE



Knowing when to replace your guard is no longer a mystery! SecurityGuard's direct-view feature lets you inspect the packing material for visual contaminants and indicates when it's time to replace the cartridge. No other guard cartridge has this convenient feature.

## EXTRA PROTECTION



SecurityGuard offers the option of stacking two cartridges in the same holder, using the simple stacking ring provided. Extra length provides extra protection. When the first cartridge becomes exhausted, contaminants are retained by the second cartridge.

## VERSATILITY



One direct-connect holder conveniently finger-tightens into virtually any brand of HPLC column worldwide. How can one holder be direct-connect and universal at the same time when end-fittings have different depths? Answer- the length of the stainless steel nib at the end of the holder automatically adjusts to the precise depth of a column's endfitting. SecurityGuard's fingertight connection will withstand pressures up to 5000 psi and it features a completely inert and biocompatible flowpath.

## ACCURACY

The cartridges can be used with virtually any matching phase of virtually any brand of column without affecting efficiency, retention time or backpressure. There are 34 different phases to choose from, including cartridges for general purpose, pharmaceutical, protein and polypeptide, aqueous size exclusion, carbohydrate and organic acid applications. SecurityGuard phases can be used with columns containing 3, 3.5, 4, 5, 10, 15  $\mu\text{m}$  or larger diameter particle sizes.

Continued on p. 26

# SECURITY GUARD ORDERING INFORMATION

## Analytical Holder Assembly Kit

Part No.	Description	Unit
<b>KJO-4282</b>	Guard Cartridge Kit	ea



*Kit includes:* 1 Cartridge Holder, 3 PEEK Ferrules, 2 Stacking Rings, 2 PEEK Fingertight Male Nuts, 2 Wrenches


## Semi-Preparative and Preparative Holder for 10.0, 21.2 and 30.0 mm ID cartridges

Part No.	Description	Unit
<b>AJO-7220</b>	Holder for 10.0 mm ID cartridges	ea
<b>AJO-8223</b>	Holder for 21.2 mm ID cartridges	ea
<b>AJO-8277</b>	Holder for 30.0mm ID cartridges	ea

*Semi-Preparative*



*Preparative*



## Cartridges - General Purpose / Pharmaceutical

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-4286</b>	C18 (ODS, Octadecyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4287</b>	C18 (ODS, Octadecyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7221</b>	C18 (ODS, Octadecyl)	1.5-10	10 x 10	3/pk
<b>AJO-7839</b>	C18 (ODS, Octadecyl)	1.5-10	15 x 21.2	ea
<b>AJO-8301</b>	C18 (ODS, Octadecyl)	1.5-10	15 x 30	ea
<b>AJO-6073</b>	C12 (Dodecyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-6074</b>	C12 (Dodecyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7275</b>	C12 (Dodecyl)	1.5-10	10 x 10	3/pk
<b>AJO-7842</b>	C12 (Dodecyl)	1.5-10	15 x 21.2	ea
<b>AJO-8304</b>	C12 (Dodecyl)	1.5-10	15 x 30	ea
<b>AJO-4289</b>	C8 (Octyl, MOS)	1.5-10	4 x 2.0	10/pk
<b>AJO-4290</b>	C8 (Octyl, MOS)	1.5-10	4 x 3.0	10/pk
<b>AJO-7222</b>	C8 (Octyl, MOS)	1.5-10	10 x 10	3/pk
<b>AJO-7840</b>	C8 (Octyl, MOS)	1.5-10	15 x 21.2	ea
<b>AJO-8302</b>	C8 (Octyl, MOS)	1.5-10	15 x 30	ea
<b>AJO-4292</b>	C5 (Pentyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4293</b>	C5 (Pentyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7372</b>	C5 (Pentyl)	1.5-10	10 x 10	3/pk
<b>AJO-4298</b>	C1 (TMS)	2-9	4 x 2.0	10/pk
<b>AJO-4299</b>	C1 (TMS)	2-9	4 x 3.0	10/pk
<b>AJO-7373</b>	C1 (TMS)	2-9	10 x 10	3/pk
<b>AJO-4347</b>	Silica	—	4 x 2.0	10/pk
<b>AJO-4348</b>	Silica	—	4 x 3.0	10/pk
<b>AJO-7223</b>	Silica	—	10 x 10	3/pk
<b>AJO-7229</b>	Silica	—	15 x 21.2	ea
<b>AJO-8312</b>	Silica	—	15 x 30	ea
<b>AJO-8328</b>	HILIC	1.5-8	4 x 2.0	10/pk
<b>AJO-8329</b>	HILIC	1.5-8	4 x 3.0	10/pk
<b>AJO-4301</b>	NH <sub>2</sub> (Amino, Aminopropyl)		4 x 2.0	10/pk
<b>AJO-4302</b>	NH <sub>2</sub> (Amino, Aminopropyl)		4 x 3.0	10/pk
<b>AJO-7364</b>	NH <sub>2</sub> (Amino, Aminopropyl)		10 x 10	3/pk
<b>AJO-8162</b>	NH <sub>2</sub> (Amino, Aminopropyl)		15 x 21.2	ea
<b>AJO-8309</b>	NH <sub>2</sub> (Amino, Aminopropyl)		15 x 30	ea

## Cartridges - General Purpose / Pharmaceutical

(Continued)

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJ0-4304</b>	CN (Cyano, Cyanopropyl)		4 x 2.0	10/pk
<b>AJ0-4305</b>	CN (Cyano, Cyanopropyl)		4 x 3.0	10/pk
<b>AJ0-7313</b>	CN (Cyano, Cyanopropyl)		10 x 10	3/pk
<b>AJ0-8220</b>	CN (Cyano, Cyanopropyl)		15 x 21.2	ea
<b>AJ0-8311</b>	CN (Cyano, Cyanopropyl)		15 x 30	ea
<b>AJ0-4350</b>	Phenyl (Phenylhexyl)	1.5-10	4 x 2.0	10/pk
<b>AJ0-4351</b>	Phenyl (Phenylhexyl)	1.5-10	4 x 3.0	10/pk
<b>AJ0-7314</b>	Phenyl (Phenylhexyl)	1.5-10	10 x 10	3/pk
<b>AJ0-7841</b>	Phenyl (Phenylhexyl)	1.5-10	15 x 21.2	ea
<b>AJ0-8303</b>	Phenyl (Phenylhexyl)	1.5-10	15 x 30	ea
<b>AJ0-8326</b>	PFP(2) (Pentafluorophenylpropyl)		4 x 2.0	10/pk
<b>AJ0-8327</b>	PFP(2) (Pentafluorophenylpropyl)		4 x 3.0	10/pk
<b>AJ0-8376</b>	PFP(2) (Pentafluorophenylpropyl)		10 x 10	3/pk
<b>AJ0-8377</b>	PFP(2) (Pentafluorophenylpropyl)		15 x 21.2	ea
<b>AJ0-8378</b>	PFP(2) (Pentafluorophenylpropyl)		15 x 30	ea
<b>AJ0-4307</b>	SCX (SA, Strong Cation Exchanger)		4 x 2.0	10/pk
<b>AJ0-4308</b>	SCX (SA, Strong Cation Exchanger)		4 x 3.0	10/pk
<b>AJ0-7369</b>	SCX (SA, Strong Cation Exchanger)		10 x 10	3/pk
<b>AJ0-4310</b>	SAX (SA, Strong Cation Exchanger)		4 x 2.0	10/pk
<b>AJ0-4311</b>	SAX (SA, Strong Cation Exchanger)		4 x 3.0	10/pk
<b>AJ0-7370</b>	SAX (SA, Strong Cation Exchanger)		10 x 10	3/pk
<b>AJ0-5808</b>	RP-1(Reversed Phase Polymer)		4 x 2.0	10/pk
<b>AJ0-5809</b>	RP-1(Reversed Phase Polymer)		4 x 3.0	10/pk
<b>AJ0-7368</b>	RP-1(Reversed Phase Polymer)		10 x 10	3/pk
<b>AJ0-8358</b>	RP-1(Reversed Phase Polymer)		15 x 21.2	ea
<b>AJ0-6075</b>	Polar-RP (Ether-linked Phenyl)		4 x 2.0	10/pk
<b>AJ0-6076</b>	Polar-RP (Ether-linked Phenyl)		4 x 3.0	10/pk
<b>AJ0-7276</b>	Polar-RP (Ether-linked Phenyl)		10 x 10	3/pk
<b>AJ0-7845</b>	Polar-RP (Ether-linked Phenyl)		15 x 21.2	ea
<b>AJ0-8307</b>	Polar-RP (Ether-linked Phenyl)		15 x 30	ea
<b>AJ0-7556</b>	Fusion-RP (C18 Polar Embedded)		4 x 2.0	10/pk
<b>AJ0-7557</b>	Fusion-RP (C18 Polar Embedded)		4 x 3.0	10/pk
<b>AJ0-7558</b>	Fusion-RP (C18 Polar Embedded)		10 x 10	3/pk
<b>AJ0-7844</b>	Fusion-RP (C18 Polar Embedded)		15 x 21.2	ea
<b>AJ0-8306</b>	Fusion-RP (C18 Polar Embedded)		15 x 30	ea
<b>AJ0-7510</b>	AQ C18 (Polar Endcapped C18)		4 x 2.0	10/pk
<b>AJ0-7511</b>	AQ C18 (Polar Endcapped C18)		4 x 3.0	10/pk
<b>AJ0-7512</b>	AQ C18 (Polar Endcapped C18)		10 x 10	3/pk
<b>AJ0-7843</b>	AQ C18 (Polar Endcapped C18)		15 x 21.2	ea
<b>AJ0-8305</b>	AQ C18 (Polar Endcapped C18)		15 x 30	ea
<b>AJ0-7596</b>	Gemini C18 (TWIN Technology)		4 x 2.0	10/pk
<b>AJ0-7597</b>	Gemini C18 (TWIN Technology)		4 x 3.0	10/pk
<b>AJ0-7598</b>	Gemini C18 (TWIN Technology)		10 x 10	3/pk
<b>AJ0-7846</b>	Gemini C18 (TWIN Technology)		15 x 21.2	ea
<b>AJ0-8308</b>	Gemini C18 (TWIN Technology)		15 x 30	ea
<b>AJ0-8367</b>	Gemini-NX (C18 TWIN-NX Technology)		4 x 2.0	10/pk
<b>AJ0-8368</b>	Gemini-NX (C18 TWIN-NX Technology)		4 x 3.0	10/pk
<b>AJ0-8369</b>	Gemini-NX (C18 TWIN-NX Technology)		10 x 10	3/pk
<b>AJ0-8370</b>	Gemini-NX (C18 TWIN-NX Technology)		15 x 21.2	ea
<b>AJ0-8371</b>	Gemini-NX (C18 TWIN-NX Technology)		15 x 30	ea
<b>AJ0-7914</b>	Gemini C6-Phenyl (TWIN Tech.)		4 x 2.0	10/pk
<b>AJ0-7915</b>	Gemini C6-Phenyl (TWIN Tech.)		4 x 3.0	10/pk
<b>AJ0-8134</b>	Oligo-RP (C18 TWIN Technology)		4 x 2.0	10/pk
<b>AJ0-8135</b>	Oligo-RP (C18 TWIN Technology)		4 x 3.0	10/pk
<b>AJ0-8136</b>	Oligo-RP (C18 TWIN Technology)		10 x 10	3/pk
<b>AJ0-8210</b>	Oligo-RP (C18 TWIN Technology)		15 x 21.2	ea
<b>AJ0-8310</b>	Oligo-RP (C18 TWIN Technology)		15 x 30	ea
<b>AJ0-8324</b>	Oligo-WAX (WA, Weak Anion Exchanger)		4 x 3.0	10/pk
<b>AJ0-8325</b>	Oligo-WAX (WA, Weak Anion Exchanger)		10 x 10	3/pk
<b>AJ0-8420</b>	Oligo-WAX (WA, Weak Anion Exchanger)		15 x 30	ea

## SECURITY GUARD ORDERING INFORMATION (CONTINUED)

### Cartridges for Protein/Polypeptide Reversed Phase

For use with all silica columns for separation of proteins and peptides, such as Jupiter® (Phenomenex); Vydac® 218TP, 214TP (Alltech Associates, Inc.); SynChropak® 300 C18, C4 (Eprogen, Inc.); Nucleosil® 300 Å C18, C4; HYPERSIL® 300 Å and all other widepore or 300 Å brands.

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-4320</b>	Widepore C18 (ODS)	1.5-10	4 x 2.0	10/pk
<b>AJO-4321</b>	Widepore C18 (ODS)	1.5-10	4 x 3.0	10/pk
<b>AJO-7224</b>	Widepore C18 (ODS)	1.5-10	10 x 10	3/pk
<b>AJO-7230</b>	Widepore C18 (ODS)	1.5-10	15 x 21.2	ea
<b>AJO-8313</b>	Widepore C18 (ODS)	1.5-10	15 x 30	ea
<b>AJO-4326</b>	Widepore C5 (Pentyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4327</b>	Widepore C5 (Pentyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7371</b>	Widepore C5 (Pentyl)	1.5-10	10 x 10	ea
<b>AJO-4329</b>	Widepore C4 (Butyl)	1.5-10	4 x 2.0	10/pk
<b>AJO-4330</b>	Widepore C4 (Butyl)	1.5-10	4 x 3.0	10/pk
<b>AJO-7225</b>	Widepore C4 (Butyl)	1.5-10	10 x 10	3/pk
<b>AJO-7231</b>	Widepore C4 (Butyl)	1.5-10	15 x 21.2	ea
<b>AJO-8314</b>	Widepore C4 (Butyl)	1.5-10	15 x 30	ea

### Cartridges for Silica GFC (Aqueous SEC)

For use with all silica GFC columns, such as BioSep™ (Phenomenex); ZORBAX® GF-series; Bio-Sil® (Bio-Rad)

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-4487</b>	GFC-2000	2-7.5	4 x 3.0	10/pk
<b>AJO-7365</b>	GFC-2000	2-7.5	10 x 10	3/pk
<b>AJO-8588</b>	GFC-2000	2-7.5	15 x 21.2	ea
<b>AJO-4488</b>	GFC-3000	2-7.5	4 x 3.0	10/pk
<b>AJO-7366</b>	GFC-3000	2-7.5	10 x 10	3/pk
<b>AJO-8589</b>	GFC-3000	2-7.5	15 x 21.2	ea
<b>AJO-4489</b>	GFC-4000	2-7.5	4 x 3.0	10/pk
<b>AJO-7367</b>	GFC-4000	2-7.5	10 x 10	3/pk
<b>AJO-8590</b>	GFC-4000	2-7.5	15 x 21.2	ea

### Cartridges for Chiral

For use with chiral columns, such as Lux™ Cellulose-1, -2, -3, -4, & Amylose-2 (Phenomenex); CHIRALCEL® OD-H®, CHIRALCEL® OJ-H®, & CHIRALPAK® AD®-H (DAICEL Chemical Industries Ltd.)

Part No.	Material Description*	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJO-8402</b>	Lux Cellulose-1	2-9	4 x 2.0	10/pk
<b>AJO-8403</b>	Lux Cellulose-1	2-9	4 x 3.0	10/pk
<b>AJO-8404</b>	Lux Cellulose-1	2-9	10 x 10	3/pk
<b>AJO-8405</b>	Lux Cellulose-1	2-9	15 x 21.2	ea
<b>AJO-8406</b>	Lux Cellulose-1	2-9	15 x 30	ea
<b>AJO-8398</b>	Lux Cellulose-2	2-9	4 x 2.0	10/pk
<b>AJO-8366</b>	Lux Cellulose-2	2-9	4 x 3.0	10/pk
<b>AJO-8399</b>	Lux Cellulose-2	2-9	10 x 10	3/pk
<b>AJO-8400</b>	Lux Cellulose-2	2-9	15 x 21.2	ea
<b>AJO-8401</b>	Lux Cellulose-2	2-9	15 x 30	ea
<b>AJO-8471</b>	Lux Amylose-2	2-9	4 x 2.0	10/pk
<b>AJO-8470</b>	Lux Amylose-2	2-9	4 x 3.0	10/pk
<b>AJO-8472</b>	Lux Amylose-2	2-9	10 x 10	3/pk

(Continued on next page)

## Cartridges for Chiral (cont'd)

Part No.	Material Description*	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJ0-8473</b>	Lux Amylose-2	2-9	15 x 21.2	ea
<b>AJ0-8474</b>	Lux Amylose-2	2-9	15 x 30	ea
<b>AJ0-8621</b>	Lux Cellulose-3	2-9	4 x 2.0	10/pk
<b>AJ0-8622</b>	Lux Cellulose-3	2-9	4 x 3.0	10/pk
<b>AJ0-8623</b>	Lux Cellulose-3	2-9	10 x 10.0	3/pk
<b>AJ0-8624</b>	Lux Cellulose-3	2-9	15 x 21.2	ea
<b>AJ0-8625</b>	Lux Cellulose-3	2-9	15 x 30.0	ea
<b>AJ0-8626</b>	Lux Cellulose-4	2-9	4 x 2.0	10/pk
<b>AJ0-8627</b>	Lux Cellulose-4	2-9	4 x 3.0	10/pk
<b>AJ0-8628</b>	Lux Cellulose-4	2-9	10 x 10.0	3/pk
<b>AJ0-8629</b>	Lux Cellulose-4	2-9	15 x 21.2	ea
<b>AJ0-8630</b>	Lux Cellulose-4	2-9	15 x 30.0	ea

\* Lux Cellulose-1 is cellulose tris(3,5-dimethylphenylcarbamate)  
 Lux Cellulose-2 is cellulose tris(3-chloro-4-methylphenylcarbamate)  
 Lux Amylose-2 is amylose tris(5-chloro-2-methylphenylcarbamate)  
 Lux Cellulose-3 is cellulose tris(4-methylbenzoate)  
 Lux Cellulose-4 is cellulose tris(4-chloro-3-methylphenylcarbamate)

## Cartridges for Carbohydrate / Organic Acid

For Organic acid and carbohydrate analysis, such as Rezex™ (Phenomenex); Aminex® (Bio-Rad); Interaction; Sugar-Pak™ (Waters).

Part No.	Material Description	pH Stability	Dimensions L x ID(mm)	Unit
<b>AJ0-4490</b>	Carbo-H <sup>+</sup>	1 - 8	4 x 3.0	10/pk
<b>AJ0-4491</b>	Carbo-Ag <sup>+</sup> *	Neutral	4 x 3.0	10/pk
<b>AJ0-4492</b>	Carbo-Pb <sup>+2</sup>	Neutral	4 x 3.0	10/pk
<b>AJ0-4493</b>	Carbo-Ca <sup>+2</sup>	Neutral	4 x 3.0	10/pk

\*For use with saccharide and oligosaccharide columns in Ag<sup>+</sup> form.

## Replacement Parts

Part No.	Description	Unit
<b>AJ0-4283</b>	PEEK Ferrules	3/pk
<b>AJ0-4285</b>	Stacking Rings	2/pk
<b>AQ0-1389</b>	PEEK Fingertight Male Nuts	10/pk
<b>AJ0-4284</b>	Security Guard Wrenches	2/pk
<b>AQ0-8374</b>	PREP Coupler, SS w/ PEEK Ferrule Inserts 10-32 Threads, 1/16 in. OD x 0.020 in. ID	ea
<b>AQ0-8375</b>	Replacement Ferrule Inserts, for PREP Coupler, PEEK, 0.020 in. ID	10/pk
<b>AQ0-8222</b>	PREP Replacement O-Rings, Kalrez® For 15 x 21.2 mm SG Holder, Size 2-021	2/pk
<b>AQ0-8318</b>	PREP Replacement O-Rings, Kalrez® For 15 x 30 mm SG Holder, Size 2-025	2/pk

# HPLC SYSTEM TEST KIT



- Diagnose hardware problems rapidly and easily
- Avoid unnecessary and costly system repairs
- Convenient benchmark testing of HPLC systems using a C18 column standard
- Test system setup and hardware connections
- Quickly isolate method development problems
- Reduce instrument down time

Each kit contains the following:

1. Phenomenex 5  $\mu$ m C18, 50 x 4.6 mm HPLC column
2. Five vials of Isocratic Test Mix
3. Five vials of Gradient Test Mix



## ORDERING INFORMATION

Part No.	Description	Unit
CHO-1684	HPLC System Test Kit, Reversed Phase, includes: C18 column, isocratic and gradient test mixes	ea
CHO-1685	Isocratic Test Mix	5/pk
CHO-1686	Gradient Test Mix	5/pk

## COLUMN PERFORMANCE CHECK STANDARDS



- Convenient way to check column performance
- Affordable and easy to use

Phenomenex offers a comprehensive line of column performance check standards to help you evaluate column performance. We recommend using the check standards to verify performance of all columns upon receiving them and periodically over the lifetime of the column. Test conditions are located in the column jacket.

### NORMAL PHASE

Part No. **ALO-3033**

(For Si, NH<sub>2</sub>, NO<sub>2</sub>, Alumina, PAC, and Luna CN)

**Unit quantity:** 2 mL

**Contains:** Meta-xylene, Nitrobenzene

### REVERSED PHASE 1

Part No. **ALO-3034**

(For C1, C18, CN and Phenyl)

**Unit quantity:** 2 mL

**Contains:** Uracil, Benzamide, Benzophenone, Biphenyl

### REVERSED PHASE 2

Part No. **ALO-3045**

(For Prodigy C8, ODS(2), ODS(3); Luna C5, C8, C18, Phenyl-Hexyl, PFP(2); Jupiter C4, C5, C18, Proteo; Columbus C8, C18; Aqua; Synergi; PhenoSphere-NEXT C8, C18; Gemini C18, C6-Phenyl; Gemini-NX C18; Clarity Oligo-RP; Kinetex C18, PFP)

**Unit quantity:** 2 mL

**Contains:** Uracil, Acetophenone, Toluene, Naphthalene

*(Please refer to the QC Test Data for specific test conditions for Jupiter and Luna)*



---

# COLUMN PERFORMANCE CHECK STANDARDS (CONTINUED)

---

## **HILIC PHASE** Part No. **ALO-8317**

(For Luna HILIC; Kinetex HILIC)

**Unit quantity:** 2 mL

**Contains:** Toluene, Uracil, Cytosine

---

## **CARBOHYDRATE MIX 1** Part No. **ALO-3035**

(For Rezex RNM, RAM and other carbohydrate analysis columns)

**Unit quantity:** 2 mL

**Contains:** Maltotriose Hydrate, Maltose, Ribitol

---

## **CARBOHYDRATE MIX 2** Part No. **ALO-3036**

(For Rezex RPM and other carbohydrate analysis columns)

**Unit quantity:** 2 mL

**Contains:** Melezitose, Glucose, Fructose, Ribitol

---

## **CARBOHYDRATE MIX 3** Part No. **ALO-3037**

(For Rezex RCM, RCU, and other carbohydrate analysis columns)

**Unit quantity:** 2 mL

**Contains:** Melezitose, Maltose, Glucose, Mannose,  
Fructose, Ribitol

---

## **OLIGOSACCHARIDE STANDARD** Part No. **ALO-3038**

(For Rezex RSO, RNO, and other oligosaccharide analysis columns)

**Unit quantity:** 2 mL

**Contains:** Light corn syrup

---

## **ORGANIC ACID STANDARD** Part No. **ALO-3039**

(For Rezex ROA and other organic acid analysis columns)

**Unit quantity:** 2 mL

**Contains:** Oxalic acid, Succinic acid, Citric acid,  
Formic acid, Tartaric acid, Acetic acid

---

## **CATION-EXCHANGE** Part No. **ALO-3040**

(For SCX, SA, CM)

**Unit quantity:** 2 mL

**Contains:** Uracil, Cytosine

---

## **ANION-EXCHANGE** Part No. **ALO-3041**

(For SAX, SB, DEAE, PEI)

**Unit quantity:** 2 mL

**Contains:** Uridine, UMP

---

# COLUMN PERFORMANCE CHECK STANDARDS (CONTINUED)

---

## AQUEOUS SEC 1

Part No. **ALO-3042**

(For BioSep-SEC-S and other protein SEC columns)

**Unit quantity:** Dry; Reconstituted to 2 mL**Contains:** Bovine thyroglobulin  
Human gamma globulin  
Ovalbumin  
Myoglobin  
Uridine*(reconstitute with 1 mL of 100 mM Sodium Phosphate pH=6.8)*

---

## AQUEOUS SEC 2

Part No. **ALO-3043**

(For PolySep GFC-P and other aqueous-soluble analysis columns)

**Unit quantity:** 2 mL**Contains:** Ethylene Glycol

---

## STAR-ION A300

Part No. **ALO-3420****Unit quantity:** 2 mL**Contains:**

<u>Conc. (mg/mL)</u>			
Fluoride	5	Nitrite	20
Nitrate	20	Sulfate	20
Chloride	10	Bromide	20
Phosphate	30		

---

## POLYMERX RP-1

Part No. **ALO-7260****Unit quantity:** 2 mL**Contains:**

<u>Conc. (mg/mL)</u>	
Cytosine	13
Uracil	13
Uridine	33

---

## ONYX MONOLITHIC REVERSED PHASE

Part No. **ALO-7836****Unit quantity:** 2 mL**Contains:**

<u>Conc. (µg/mL)</u>	
Thiourea	10
Progesterone	100
Anthracene	10

---

## ONYX MONOLITHIC NORMAL PHASE

Part No. **ALO-7835****Unit quantity:** 2 mL**Contains:**

<u>Conc. (µg/mL)</u>	
Toluene	21.75
Nitrobenzene	150
2-Nitroanisol	0.18

---

## COLUMN PERFORMANCE CHECK STANDARDS (CONTINUED)

### **CHIRAL TEST MIX NO. 1**

Part No. **ALO-3046**

Applicable to the following Chirex columns:  
3001, 3005

**Unit quantity:** 2 mL

**Contains:** 1. S-(+)-2,2,2-trifluoro-1-(9-anthryl)  
ethanol CAS [60646-30-2]  
2. R-(-)-2,2,2-trifluoro-1-(9-anthryl)  
ethanol CAS [53531-34-3]

### **CHIRAL TEST MIX NO. 2**

Part No. **ALO-3047**

Applicable to the following Chirex columns:  
3010, 3011, 3012

**Unit quantity:** 2 mL

**Contains:** N-dansyl-DL-valine  
(cyclohexylammonium salt)  
CAS[84540-67-0]

### **CHIRAL TEST MIX NO. 3**

Part No. **ALO-3048**

Applicable to the following Chirex columns:  
3014, 3017, 3018, 3019, 3020, 3022

**Unit quantity:** 2 mL

**Contains:** 1.(R)-(-)-N-(3,5-Dinitrobenzoyl)- $\alpha$ -  
methylbenzylamine CAS [69632-32-2]  
2.(S)-(-)-N-(3,5-Dinitrobenzoyl)- $\alpha$ -  
methylbenzylamine CAS[69632-31-1]

### **CHIRAL TEST MIX NO. 4**

Part No. **ALO-3049**

Applicable to the following Chirex column:  
3126

**Unit quantity:** 2 mL

**Contains:** DL-Aspartic Acid CAS [617-45-8]

### **CHIRAL TEST MIX NO. 5**

Part No. **ALO-8412**

Applicable to the following Lux columns:  
Lux Cellulose -1,-2,-3,-4, Lux Amylose-2

**Unit quantity:** 2 mL

**Contains:** Trans-Stilbene oxide CAS [1439-07-2]

## PART XII - SOLID PHASE EXTRACTION (SPE)

Increase column and instrument life by injecting samples cleaned-up with Strata®.

### STRATA™-X Polymeric Sorbents

Tubes and 96-Well Plates

- Deconditioning Resistant
- Low Elution Volumes
- High Analyte Capacity



**Strata™-X** for simplified cleanup of polar and non-polar compounds

**Strata™-X-C** for selective extraction of basic compounds

**Strata™-X-CW** for bases (including quaternary amines)

**Strata™-X-A** for cleanup of weak acids

**Strata™-X-AW** for acids

### STRATA® Traditional Sorbents

Tubes and 96-Well Plates

- Optimal Flow
- Lot-to-Lot Reproducibility
- Wide Range of Selectivity
- Available chemistries include:  
C18-E, C18-U, C18-T, C8,  
Phenyl, SDB-L, CN, Si-1, WCX, FI-PR, NH<sub>2</sub>, SAX, SCX,  
Melamine



### STRATA® Flash Sorbents

- Polar & Non-polar Phases
- Narrow Particle Range Distribution
- Can be used for Direct Scale-up



**Strata® Giga™ Tubes** available in 12, 20, 60 & 150 mL formats

**Sepra™ Bulk** available in gram to multi-kilogram quantities

### STRATA® On-line Cartridges

- Rapid Extraction and Concentration
- Direct Inject Analysis
- Easily Automated



**Strata™-X** for polar and non-polar compounds

**Strata™-X-C** for weak bases

**Strata™-X-CW** for strong bases

**Strata® C18** for non-polar compounds

**Strata® C8** for compounds of intermediate polarity

Aqua, Gemini, Jupiter, Luna, and Strata are registered trademarks of Phenomenex, Inc. Axia, BioSep, Columbus, Giga, Kinetex, Lux, Onyx, Phenex, Phenogel, PhenoSphere, PolymerX, PolySep, Prodigy, Rezex, SecurityGuard, Sepra, Star-Ion, Strata-X, Synergi, ThermoSphere, TWIN-NX and TWIN Technology are trademarks of Phenomenex, Inc. Alltima is a trademark of Alltech Associates, Inc. Bio-Sil and Aminex are registered trademarks of Bio-Rad Laboratories. Chirex is a registered trademark of Chirex, Inc. licensed to Phenomenex. Cyrolite is a trademark of CY/RO Industries. HYPERSIL is a registered trademark of Thermo Hypersil-Keystone. Nucleosil is a registered trademark of Macherey-Nagel. Sugar-Pak is a trademark of Waters Corporation. Symmetry is a registered trademark of Waters Corporation. SynChropak is a registered trademark of Eprogen, Inc. Teflon is a registered trademark of E.I. du Pont de Nemours and Co. Vydac is a registered trademark of Alltech Associates, Inc. ZORBAX is a registered trademark of Agilent Technologies. © 2010 Phenomenex. All rights reserved. Onyx is a product based on monolithic technology under license from Merck KGaA, Darmstadt, Germany.

**Subject to Phenomenex Standard Terms & Conditions which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).**

## PART XIII - HPLC ACCESSORIES

### ACCESSORIES

- Backpressure Regulators
- Biocompatible / Metal-free products
- Connectors and Splitters
- Filtration Products
- Injectors and Injector Loops
- Membrane Filters
- Mobile Phase Handling Devices
- Polymer Calibration Standards / Kits
- Rotor Seals, Stators, etc.
- Solvent Reservoir and Reagent Bottles
- SPE Consumables, Tube & Plate Manifolds
- Switching Valves
- Syringes
- Syringe Filters
- Tools
- Tubing, Fittings, Frits and Unions
- Valves (Injection, Switching)
- Vials, Caps and Septa

---


### EQUIPMENT

- Column Chiller-Heater
- Column Heater
- Column Selectors
- Degasser
- Fluid Processors
- Mobile Phase Recycler
- Temperature Controllers

*For ordering and additional information, please contact your Phenomenex Technical Consultant.*

### SINGLE COLUMN HEATER THERMASPHERE™ TS-130

- Compact, low-cost heater precisely controls temperature from 25-90 °C
- Improves reproducibility and chromatographic results
- Reduces analyte identification errors
- Improves baseline and overall detector performance
- Improves peak efficiency and analyte quantitation (especially at low levels)
- Improves the ruggedness of the separation (within-lab and lab-to-lab)



---

### ORDERING INFORMATION

#### ThermaSphere™ TS-130

Part No.	Description
<b>EHO-7057</b>	ThermaSphere TS-130 HPLC Column Heater 25-90 °C, 95 to 265 VAC, 50/60 Hz
<b>EHO-7058</b>	Stand for ThermaSphere TS-130 HPLC Column Heater

*More Accessories available. See Phenomenex Catalog for details.*

# PHENOMENEX WARRANTY

Phenomenex products are warranted to meet the stated performance and quality and to be free of defects in material and workmanship. They are not warranted, nor does Phenomenex assume liability, if misused. **NO OTHER WARRANTY OR REPRESENTATION IS IMPLIED OR EXPRESSED BY PHENOMENEX FOR ITS PRODUCTS WITH RESPECT TO MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR ANY OTHER MATTER. PHENOMENEX SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, OR COMPENSATORY DAMAGE ARISING FROM THE USE OF, OR IN CONJUNCTION WITH, ITS PRODUCTS.** The maximum liability which can be assumed by Phenomenex for breach of warranty shall be the invoice price of the product.

---

## SPECIFIC WARRANTIES ON HPLC COLUMNS

Phenomenex warrants its quality columns in accordance with the following terms and conditions. Phenomenex will repack, replace, or refund charges on any column (at our discretion), at no cost if a column fails to perform satisfactorily. Columns being returned must have prior return authorization granted by Phenomenex. Defective products must be accompanied by a written explanation of failure. Approval is subject to the following exclusions:

- All columns must be tested upon receipt and all deficiencies must be reported to Phenomenex no later than 15 days after the date of receipt of the column.
- Maximum warranty period is limited to 90 days on HPLC columns unless previously agreed upon. However, **COLUMNS MAY NOT BE RETURNED FOR REFUND OR CREDIT AFTER 45 DAYS AND WITHOUT PRIOR AUTHORIZATION.**
- Removal of column end-fittings automatically voids column warranty.
- Column performance warranty is limited to the conditions of the original test chromatograms.
- Physical damage to the column due to misuse, abuse, or mishap, including mechanical shock.
- Chemical damage to the packing material due to operation at incorrect chemical conditions, temperatures, or pressures.
- Failure due to high backpressures caused by improper solvent or sample filtration practices causing particulate build-up or precipitation in the column or end-fitting.
- Incorrect selection of packing material made by customer for their particular use or incompatibility of equipment, etc.
- For products supplied by but not manufactured by Phenomenex, the warranty is limited by the terms of the original manufacturer's warranty.





**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country,  
contact Phenomenex USA, International Department by telephone, fax or email: [international@phenomenex.com](mailto:international@phenomenex.com).



*...breaking with tradition<sup>sm</sup>*

[www.phenomenex.com](http://www.phenomenex.com)

[support@phenomenex.com](mailto:support@phenomenex.com)