

## A Simplified and Streamlined Approach to Solid Phase Extraction using Strata™-X-Drug B for SAMHSA's Low Cutoff Drugs of Abuse Panel

Michael Rummel, Matthew Trass, Seyed Sadjadi, and Erica Pike  
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Utilizing a single Solid Phase Extraction (SPE) sorbent and only 3 extraction methods, 11 drugs from 6 different Substance Abuse and Mental Health Services Administration (SAMHSA) regulated drug classes were effectively extracted and analyzed by LC/MS and GC/MS. The SPE sorbent, Strata-X-Drug B, was specifically designed and QC tested for the forensic toxicology environment and as a result was able to provide consistent high recoveries of each drug analyzed that conformed to the new lower cutoff levels that SAMHSA has recently enforced.

### Introduction

Workplace drug testing has become routine in the job application process. Corporations are also adopting workplace substance abuse policies which require employees to submit to random drug testing. Drug testing is not only increasing in the workplace but it also plays a large role in rehabilitation facilities and law enforcement agencies.<sup>1</sup> With drug testing on the rise, it has become difficult to establish analytical methods for each of the many drugs that are regulated by SAMHSA because each drug class is vastly different from the next and the regulated cutoff levels are getting lower and lower (**Table 1**). The goal of this work was to develop extraction methods using a single sorbent for all SAMHSA regulated drugs that meet the most recent SAMHSA dictated cutoff levels. Not only was it important to use a single sorbent but it was also important to develop as few extraction methods as possible in order to reduce user error and technician training times while still producing consistent, optimal results. These requirements led to the creation of a Solid Phase Extraction (SPE) sorbent, Strata-X-Drug B, which is able to extract all drugs regulated by SAMHSA in 3 short methods. To ensure consistent recoveries and reliable results, Strata-X-Drug B is also specifically quality control tested for the forensic toxicology environment by testing recoveries with drugs of abuse probes from urine as well as testing the sorbent to be certain that it does not promote interconversion of norcodeine and normorphine to their parent compounds, which can sometimes happen with strong cation-exchange SPE sorbents.

**Table 1.**  
SAMHSA Regulated Cutoff Levels<sup>2</sup>

Class	Analyte	Cutoff
Marijuana Metabolites	Delta-9-tetrahydrocannabinol-9-carboxylic acid (THCA)	15 ng/mL
Cocaine Metabolites	Benzoylcegonine	100 ng/mL
Opiate Metabolites	Codeine	2,000 ng/mL
	Morphine	2,000 ng/mL
6-MAM	6-Acetylmorphine	10 ng/mL
PCP	Phencyclidine	25 ng/mL
Amphetamines	Amphetamine	250 ng/mL
	Methamphetamine	250 ng/mL
MDMA	Methylenedioxyamphetamine (MDMA)	250 ng/mL
	Methylenedioxyamphetamine (MDA)	250 ng/mL
	Methylenedioxyethylamphetamine (MDEA)	250 ng/mL

### Experimental Conditions

Several drug classes under the expanded SAMHSA regulations were analyzed by LC/MS and GC/MS. The drug classes include opiates, 6-MAM, PCP, amphetamines, marijuana metabolites, and cocaine metabolites. Standards were obtained from Cerilliant® Corporation (Round Rock, TX). Urine samples were prepared for each drug class and were spiked at 40, 100, and 125 % of the SAMHSA cutoff level. The prepared urine samples were then subjected to various pre-treatment steps, depending on the drug that was spiked into each (**Table 2**). Pre-treated samples were then subjected to SPE using 60 mg/6 mL Strata-X-Drug B tubes (Phenomenex, Inc. Torrance, CA, USA) as specified in **Table 3**. After extraction, each of the 6 drug classes were analyzed by either GC/MS or LC/MS under the conditions specified in **Figures 1-6**. GC/MS analysis was performed on a Zebtron™ ZB-Drug-1 column (Phenomenex, Inc. Torrance, CA, USA) with the MS operated in SIM mode while LC/MS analysis was performed either using a Kinetex® C18 or a PFP core-shell HPLC column (Phenomenex, Inc. Torrance, CA, USA) with the MS operating in ESI+ mode.

**Table 2.**  
Sample Pre-treatments

Opiates	To each 2 mL urine sample, add 500 µL of concentrated Hydrochloric acid. Heat at 90 °C for 2 hours. Add 2 mL of 200 mM Sodium acetate buffer (pH 4.0). Add 1 mL of 6 N Potassium hydroxide then vortex. Centrifuge for 5 minutes at 5000 rpm. Verify that pH is between 4.0-6.0.
6-MAM	To each 2 mL urine sample, add 1000 µL of β-glucuronidase solution (contains 5000 F units/mL <i>Patella vulgata</i> in 100 mM Acetate buffer, pH 5.0) and vortex. Hydrolyze for 3 hours at 60 °C. Let cool and add 1000 µL of 100 mM Phosphate buffer (pH 6.0). Verify that pH is between 5.5 and 6.5. Centrifuge for 5 minutes at 5000 rpm and discard pellet.
PCP	To each 2 mL urine sample, add 2 mL of 100 mM Sodium acetate buffer (pH 5.0) then vortex. Verify that pH is between 4.0 and 6.0.
Amphetamines	To each 2 mL urine sample, add 1000 µL of 100 mM Phosphate buffer (pH 6.0) and 1000 µL of 0.35 M Sodium periodate. Vortex and incubate at room temperature for 25 minutes. Verify that the pH is between 5.5 and 6.5.
Marijuana Metabolites	To each 2 mL urine sample, add 100 µL of 11.8 N Potassium hydroxide. Vortex then incubate at 60 °C for 20 minutes. Cool and add ~450 µL glacial acetic acid then vortex. Verify that the pH is between 4.0 and 6.0.
Cocaine Metabolites	To each 2 mL urine sample, add 2 mL of 100 mM Sodium acetate buffer (pH 5.0) then vortex. Verify that the pH is between 4.0 and 6.0.

# TN-0041 APPLICATIONS

**Table 3.**

Protocols


Strata™-X-Drug B 60 mg/6 mL (part no. 8B-S128-UCH)

	Opiates, 6-MAM, PCP, and Amphetamines	Marijuana Metabolites	Cocaine Metabolites
Condition/Equilibrate	No Conditioning Required		
Load	Pre-treated urine sample	Pre-treated urine sample	Pre-treated urine sample
Wash 1	2 mL of 100 mM Sodium acetate buffer (pH 5.0)	2 mL of 100 mM Sodium acetate (pH 5.0)	2 mL of 0.1 N Hydrochloric acid
Wash 2	2 mL Methanol	2 mL of (30:70) Acetonitrile/100 mM Sodium acetate (pH 5.0)	2 mL Methanol
Dry	10 minutes under full vacuum	15 minutes under full vacuum	10 minutes under full vacuum
Elute	2 mL of Ethyl acetate/ Isopropanol/ Ammonium hydroxide (70:20:10)	2 mL of Ethyl acetate/ Isopropanol (85:15)	2 mL of Ethyl acetate/ Isopropanol/ Ammonium hydroxide (70:20:10)

**Figure 1.**

GC Analysis of Opiates


Extracted samples were evaporated to dryness at 50 °C then reconstituted in 50 µL of propionic anhydride (PIA) and 50 µL of acetonitrile. Samples were heated for 30 minutes at 60 °C and were then evaporated to dryness at 50 °C. Samples were again reconstituted in 100 µL of ethyl acetate.



**Figure 2.**

HPLC Analysis of 6-MAM


Extracted samples were evaporated to dryness at 50 °C then reconstituted in 1 mL of 10 % Methanol / 90 % Formic acid (0.1 %).



**Figure 3.**

HPLC Analysis of PCP


Extracted samples were evaporated to dryness at 50 °C then reconstituted in 1 mL of 25 % Methanol / 75 % Formic acid (0.1 %).



**Figure 4.**

HPLC Analysis of Amphetamines


300 µL of 0.5 N Methanolic hydrochloride (0.5 N HCl in Methanol) was added to each sample. Samples were then evaporated to dryness at < 35 °C then reconstituted with 1 mL of 10 % Methanol / 90 % Formic acid (0.1 %). Before injecting onto the LC/MS, samples were diluted by a factor of 20 to bring the concentration into a suitable range for analysis.



# TN-0041 APPLICATIONS


**Figure 5.**  
GC Analysis of Marijuana Metabolites

Extracted samples were evaporated to dryness at 50 °C then reconstituted in 25 µL of ethyl acetate and 50 µL of BSTFA. Samples were heated for 20 minutes at 70 °C.



**Figure 6.**  
HPLC Analysis of Cocaine Metabolites

Extracted samples were evaporated to dryness at 50 °C then reconstituted in 1 mL of 95 % Methanol / 5 % Formic acid (0.1 %).



## Results and Discussion

While developing the extraction methods, it was discovered that conditioning and equilibrating the Strata-X-Drug B sorbent was not necessary. The elimination of these steps reduced solvent consumption and time without compromising results. The unconditioned sorbent proved to be effective in the extraction of all 11 drugs regulated under SAMHSA by providing high recoveries and consistent results (**Table 4**). Coupled with efficient GC/MS and LC/MS analysis using either a Zebron ZB-Drug 1 GC column or a Kinetex core-shell HPLC column, it was determined that the extraction and analysis techniques were able to produce clean samples, excellent separation and quantitation, and quick analysis times which conform with SAMHSA regulations.

**Table 4.**  
Relative Recovery, RSD, and Linearity

Analyte	Relative Recovery (%)	RSD % (N=3)	Linearity
Codeine	99.86	0.43	0.999
Morphine	99.80	0.59	0.999
6-MAM	100.5	0.63	0.999
Phencyclidine	99.33	0.63	0.999
Amphetamine	101.8	1.59	0.997
Methamphetamine	97.76	1.99	0.998
MDA	96.65	3.03	0.995
MDMA	100.1	1.20	0.998
MDEA	102.2	1.93	0.999
Carboxy-THC	99.45	0.95	0.999
Benzoylecgonine	99.45	0.95	0.999

## Conclusion

Strata-X-Drug B effectively decreased extraction time by eliminating the condition and equilibration steps, reducing the extraction method steps to a load, two washes, and an elution. Strata-X-Drug B also allowed for a streamlined extraction approach by minimizing the number of different extraction protocols required for the expanded SAMHSA analysis to only 3 methods with few variations of solvents between methods. This approach drastically reduced processing time and user error allowing for increased throughput and higher quality data. It was also shown that subsequent analysis by GC/MS and LC/MS proved that both analysis techniques were acceptable under the same SPE extraction methodology. However, LC/MS proved to be ideal as it resulted in much shorter analysis times which allowed for a significant increase in sample throughput.

## References

1. [www.wikipedia.org/wiki/Drug-Testing](http://www.wikipedia.org/wiki/Drug-Testing)
2. [www.samhsa.gov](http://www.samhsa.gov)

# TN-0041 APPLICATIONS

## Australia

t: 02-9428-6444  
f: 02-9428-6445  
auinfo@phenomenex.com

## Austria

t: 01-319-1301  
f: 01-319-1300  
anfrage@phenomenex.com

## Belgium

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
beinfo@phenomenex.com

## Canada

t: (800) 543-3681  
f: (310) 328-7768  
info@phenomenex.com

## Denmark

t: 4824 8048  
f: 4810 6265  
nordicinfo@phenomenex.com

## Finland

t: (09)4789 0063  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## France

t: 01 30 09 21 10  
f: 01 30 09 21 11  
franceinfo@phenomenex.com

## Germany

t: 06021-58830-0  
f: 06021-58830-11  
anfrage@phenomenex.com

## Ireland

t: 01 247 5405  
f: +44 1625-501796  
eireinfo@phenomenex.com

## Italy

t: 051 6327511  
f: 051 6327555  
italiainfo@phenomenex.com

## Luxembourg

t: +31 (0)30-2418700  
f: +31 (0)30-2383749  
nlinfo@phenomenex.com

## Mexico

t: (55) 5018 3791  
f: (310) 328-7768  
tecnicomx@phenomenex.com

## The Netherlands

t: 030-2418700  
f: 030-2383749  
nlinfo@phenomenex.com

## New Zealand

t: 09-4780951  
f: 09-4780952  
nzinfo@phenomenex.com

## Norway

t: 810 02 005  
f: +45 4810 6265  
nordicinfo@phenomenex.com

## Puerto Rico

t: (800) 541-HPLC  
f: (310) 328-7768  
info@phenomenex.com

## United Kingdom

t: 01625-501367  
f: 01625-501796  
ukinfo@phenomenex.com

**All other countries:**  **Corporate Office USA**

t: (310) 212-0555  
f: (310) 328-7768  
info@phenomenex.com

## www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

## Ordering Information

### Strata™ -X-Drug B SPE

Sorbent Mass	Part No.	Unit
<b>Tube</b>		
10 mg	8B-S128-AAK	1 mL (100/box)
10 mg	8L-S128-AAK	1 mL (100/box)
30 mg	8B-S128-TAK	1 mL (100/box)
30 mg	8L-S128-TAK	1 mL (100/box)
30 mg	8B-S128-TBJ	3 mL (50/box)
60 mg	8B-S128-UBJ	3 mL (50/box)
60 mg	8B-S128-UCH	6 mL (30/box)
60 mg	8B-S128-UCL	6 mL (200/box)
<b>Giga™ Tube</b>		
100 mg	8B-S128-EDG	12 mL (20/box)
<b>96-Well Plate</b>		
10 mg	8E-S128-AGB	2 Plates/box
30 mg	8E-S128-TGB	2 Plates/box
60 mg	8E-S128-UGB	2 Plates/box

\* Tab-less tube

### Kinetex® Core-Shell HPLC/UHPLC Columns

#### 1.7 µm Minibore Columns (mm)

	50 x 2.1	100 x 2.1	150 x 2.1
<b>C18</b>	00B-4475-AN	00D-4475-AN	00F-4475-AN
<b>PFP</b>	00B-4476-AN	00D-4476-AN	00F-4476-AN

#### 2.6 µm Minibore Columns (mm)

	50 x 2.1	100 x 2.1	150 x 2.1
<b>C18</b>	00B-4462-AN	00D-4462-AN	00F-4462-AN
<b>PFP</b>	00B-4477-AN	00D-4477-AN	00F-4477-AN

#### 2.6 µm Solvent Saver MidBore™ Columns (mm)

	50 x 3.0	100 x 3.0	150 x 3.0
<b>C18</b>	00B-4462-Y0	00D-4462-Y0	00F-4462-Y0
<b>PFP</b>	00B-4477-Y0	00D-4477-Y0	00F-4477-Y0

#### 2.6 µm Analytical Columns (mm)

	50 x 4.6	100 x 4.6	150 x 4.6
<b>C18</b>	00B-4462-E0	00D-4462-E0	00F-4462-E0
<b>PFP</b>	00B-4477-E0	00D-4477-E0	00F-4477-E0

More dimensions and phases available, please inquire.

### Zebtron™ ZB-Drug-1 GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.18	0.18	40 to 320/340 °C	7CD-G023-08
<b>15-Meter</b>			
0.25	0.25	40 to 320/340 °C	7EG-G023-11
<b>30-Meter</b>			
0.25	0.25	40 to 320/340 °C	7HG-G023-11



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to another product of the same phase and comparable dimensions, return the product with comparative data within 45 days for a FULL REFUND.

#### Terms and Conditions

Subject to Phenomenex Standard Terms and Conditions which may be viewed at <http://www.phenomenex.com/TermsAndConditions>.

#### Trademarks

Strata-X, Giga, MidBore, and Zebtron, are trademarks of Phenomenex, Inc. Kinetex is a registered trademark of Phenomenex in the United States, European Union, and other jurisdictions. Cerilliant is a registered trademark of Cerilliant Corporation.

#### Disclaimer

Phenomenex, Inc. is in no way affiliated with Cerilliant Corporation. Comparative separations may not be representative of all applications.

Strata-X is patented by Phenomenex. U.S. Patent No. 7,119,145

© 2010 Phenomenex, Inc. All rights reserved.