TN-1151 APPLICATIONS



Simple and Fast Quantitation of Nicotinic Acid and Nicotinamide in Human Plasma by Applying Impact[™] Protein Precipitation Plate Technology with Gemini[®] 3 µm C18 HPLC Columns

Shuguang Li and Erica Pike Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Nicotinic acid and nicotinamide were extracted from human plasma by performing a rapid protein precipitation using Impact Protein Precipitation Plates followed by HPLC analysis using a Gemini 3 μ m C18 100 x 4.6 mm HPLC column and positive polarity ESI LC/MS/MS system. Impact technology offers easy, fast protein removal while providing maximized recovery of the target analytes. The Gemini 3 μ m C18 HPLC column produced excellent chromatographic resolution, sensitivity, and high peak capacities.

Introduction

Niacin (nicotinic acid) is a water-soluble vitamin that is also referred to as vitamin B3. Nicotinamide (nicotinic acid amide) is the derivative of niacin that is incorporated into the coenzymes nicotinamide adenine dinucleotide (NAD) and nicotinamide adenine dinucleotide phosphate (NADP). The nicotinamide moiety of NAD and NADP serves as an electron acceptor or donor in biological oxidation-reduction reactions catalyzed by several hundred different enzymes. Both nicotinic acid and nicotinamide are absorbed from the normal diet. Nicotinamide is the form of vitamin B3 that is commonly found in nutritional supplements and used to fortify foods. Nicotinic acid is available both over the counter and with a prescription as a cholesterol-lowering agent. Niacin deficiency can result from inadequate dietary intake of niacin and/ or the amino acid tryptophan. Niacin deficiency can affect the skin, digestive system, and the nervous system. Severe niacin deficiency is referred to as pellagra.

Materials and Methods

Sample Preparation:

Protein precipitation was performed using an Impact Protein Precipitation Plate.

Step	
1.	Place the Impact plate onto a suitable 96-well sample manifold
2.	Dispense 300 μL acetonitrile into each well of the Impact plate
3.	Add 100 μL of plasma/serum samples to each well of the Impact plate
4.	Mix 3 times by aspirating with a pipette tip
5.	Apply vacuum to filter the sample and collect the purified filtrate

After filtrate is collected, the collection plate containing the purified samples should be covered using a sealing mat. The sample is now ready to be injected onto the LC/MS/MS. If the sample will not be injected onto the LC/MS/MS immediately, transfer the filtrate to amber Verex[™] autosampler vials (ambient) to protect from light.

HPLC Conditions: Column: Gemini 3 µm C18 Dimensions: 100 x 4.6 mm Part No.: 00D-4439-E0 Mobile Phase: A: 0.1 % Formic acid in water B: Methanol Flow Rate: 600 µL/min Gradient: Time (min) % B 10 0 2.5 90 2.6 10 Δ 10 Detection: API 4000[™] MS/MS, ESI Positive (ESI+) Temperature: Ambient Injection: 2 µL purified plasma

MS/MS Conditions:

An AB SCIEX API 4000[™] triple-quadrupole tandem mass spectrometer is used for analysis equipped with an ESI probe operating in positive polarity mode. Under an MRM mode, two channels were monitored for nicotinamide and nicotinic acid (**Table 1**).

Table 1.

MRM Transitions

Peak Name	MRM Channel
Nicotinamide	123.006 ⇒ 80.100
Nicotinic acid	123.981 ⇒ 80.100

Results and Discussion


When developing a method for the analysis of nicotinamide and nicotinic acid, it was important that the method be rapid, sensitive, and accurate in order to accommodate high-throughput labs that analyze 100's to 1000's of samples each week. Traditionally, a protein precipitation step is used for fast cleanup of plasma samples. Protein precipitation is normally performed using a centrifuge tube or a 96-well collection plate, however this process requires that supernatant be collected while being careful not to disrupt pelleted protein in the bottom of the tube or collection plate. This step was greatly simplified by using Impact Protein Precipitation Plates. The Impact plate allows for the analysis of 96 samples at once, eliminates the transfer steps that are com-

TN-1151 APPLICATIONS

monly associated with protein precipitation, and can also be automated. Protein precipitation was performed within the wells of the Impact[™] plate and sample was not allowed to pass through the filter of the plate until vacuum was applied. This ensured that the precipitated protein was left within the wells of the Impact plate while protein free sample was allowed to pass through the filter and into a collection plate (**Figure 1**).


Figure 1.

Protein Precipitation Using Impact Protein Precipitation Plates




After the plasma samples were cleaned up, they were analyzed by LC/MS/MS using a Gemini® 3 µm C18 HPLC column coupled to an API 4000[™] triple-quadrapole tandem mass spectrometer. The Gemini 3 µm C18 HPLC column was chosen because it contains a unique silica-organic layer that is grafted onto the base silica which mechanically strengthens the particle while providing excellent efficiencies. Efficiency and resolution were necessary in this analysis because the chemical properties of nicotinamide and nicotinic acid are similar in that they share the same backbone, however nicotinamide contains an amide moiety while nicotinic acid has a carboxylic acid moiety (Figure 2). It was also important to minimize tailing and improve peak shape of both target compounds. Using a 100 x 4.6 mm column reduced tailing and provided the desired peak shape of each target compound. Shorter columns, such as a 50 x 4.6 mm, resulted in tailing and undesirable peak shapes indicating that a longer column is necessary to accurately separate nicotinamide and nicotinic acid.


Figure 2. Chemical Structures of Nicotinamide and Nicotinic Acid












TN-1151 APPLICATIONS

Figure 5.

Representative Standard Curve of Nicotinic Acid at a Concentration Range of 2 to 2000 ng/mL



The reproducibility of our analysis was determined by producing a standard curve of nicotinic acid at a concentration range of 2 to 2000 ng/mL. With $R^2 = 0.9996$, our method proved to be reproducible even at low levels of detection. The signal-to-noise ratio was also studied at LOD and LOQ. With a signal-to-noise ratio of 4.3 and 3.3, it was determined that we were able to reliably analyze nicotinamide and nicotinic acid from plasma/serum samples at a LOD of 2 ng/mL. The LOQ was subsequently determined to be accurate at 10 ng/mL with signal-to-noise ratios of 12.4 and 10.1.

Table 2.

Signal-to-Noise Ratio of Nicotinamide and Nicotinic Acid at LOD and LOQ

Analyte	L	DD	LOQ		
	ng/mL S/N Ratio		ng/mL	S/N Ratio	
Nicotinamide	2	4.3	10	12.4	
Nicotinic acid	2	3.3	10	10.1	

Conclusion

As the study of vitamins becomes an increasingly important factor in clinical research, high-throughput laboratories must adopt rapid and robust methods to accurately analyze and quantitate vitamins and their derivatives. With these goals in mind, we developed a method that can be easily automated, can be used in high-throughput labs, and provides both sensitivity and speed. The sample preparation step using Impact[™] Protein Precipitation Plates is simple, requires no method development, and can process 96 samples at once. The downstream LC/MS/MS analysis using a Gemini[®] 3 µm C18 HPLC column provides resolution between nicotinic acid and its derivative, nicotinamide, in under 3 minutes. The analysis is also sensitive, with an LOD at 2 ng/mL and a LOQ at 10 ng/mL for both analytes.

TN-1151 LICATION



Gemini [®] HPLC Columns										
3 μm Microbore, Minibore and Narrow Bore Columns (mm) SecurityGuard™										
Phases	50 x 1.0	20 x 2.0	30 x 2.0	50 x 2.0	100 x 2.0	150 x 2.0	50 x 3.0	100 x 3.0	150 x 3.0	4 x 2.0*
										10/pk
C18	00B-4439-A0	00M-4439-B0	00A-4439-B0	00B-4439-B0	00D-4439-B0	00F-4439-B0	00B-4439-Y0	00D-4439-Y0	00F-4439-Y0	AJ0-7596
										for ID: 2.0-3.0 mm

3 µm Analy	Guard Cartridges (mm)						
Phases	20 x 4.0	30 x 4.6	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
							10/pk
C18	00M-4439-D0	00A-4439-E0	00B-4439-E0	00D-4439-E0	00F-4439-E0	00G-4439-E0	AJ0-7597
							for ID: 3.2-8.0 mm

nzinfo@phenomenex.com

nordicinfo@phenomenex.com

nordicinfo@phenomenex.com

ukinfo@phenomenex.com

info@phenomenex.com

Corporate Office USA

info@phenomenex.com

All other countries:

New Zealand

t: 09-4780951

Norway

t: 810 02 005

Sweden t: 08 611 6950

f: +45 4810 6265

f: 01625-501796

United States t: (310) 212-0555

(310) 328-7768

t: (310) 212-0555

(310) 328-7768

United Kingdom t: 01625-501367

f: +45 4810 6265

Puerto Rico

(800) 541-HPLC

(310) 328-7768

info@phenomenex.com

09-4780952

f:

t:

f:

f:

f:

5 µm Minibore and Narrow Bore Columns (mm) SecurityGu									
Phases	30 x 2.0	50 x 2.0	150 x 2.0	250 x 2.0	50 x 3.0	100 x 3.0	150 x 3.0	250 x 3.0	4 x 2.0*
									10/pk
C18	00A-4435-B0	00B-4435-B0	00F-4435-B0	00G-4435-B0	00B-4435-Y0	00D-4435-Y0	00F-4435-Y0	00G-4435-Y0	AJ0-7596
									for ID: 2.0-3.0 mm

Luxembourg	5 µm Ana	alytical Columns	(mm)			SecurityGua	ard Cartridges (mm)
t: +31 (0)30-2418700	Phases	30 x 4.6	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	4 x 3.0*
f: +31 (0)30-2383749							10/pk
nlinfo@phenomenex.com	C18	00A-4435-E0	00B-4435-E0	00D-4435-E0	00F-4435-E0	00G-4435-E0	AJ0-7597
Mexico							for ID: 3.2-8.0 mm
t: 001-800-844-5226 f: 001-310-328-7768 tecnicomx@phenomenex.com	* SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282						
The Netherlands t: 030-2418700 f: 030-2383749 nlinfo@phenomenex.com		•				Unit	
	t: +31 (0)30-2418700 f: +31 (0)30-2383749 nlinfo@phenomenex.com Mexico t: 001-800-844-5226 f: 001-310-328-7768 tecnicomx@phenomenex.com The Netherlands t: 030-2418700	t: +31 (0)30-2418700 f: +31 (0)30-2383749 nlinfo@phenomenex.com Mexico t: 001-800-844-5226 f: 001-310-328-7768 tecnicomx@phenomenex.com The Netherlands t: 030-2418700 f: 030-2383749 plinfo@phenomenex.com	t: +31 (0)30-2418700 f: +31 (0)30-2383749 nlinfo@phenomenex.com Mexico t: 001-800-844-5226 f: 001-310-328-7768 tecnicomx@phenomenex.com The Netherlands t: 030-2418700 f: 030-2418700 f: 030-2383749 nlinfo@phenomenex.com	t: +31 (0)30-2418700 Phases 30 x 4.6 50 x 4.6 f: +31 (0)30-2383749 nlinfo@phenomenex.com Mexico C18 00A-4435-E0 00B-4435-E0 t: 001-800-844-5226 C18 00A-4435-E0 00B-4435-E0 f: 001-310-328-7768 * SecurityGuard Analytical cartridges require holder, tecnicomx@phenomenex.com * SecurityGuard Analytical cartridges require holder, t: 030-2418700 Ordering Information f: 030-2383749 Impact [™] Precipitation Products	t: +31 (0)30-2418700 Phases 30 x 4.6 50 x 4.6 100 x 4.6 f: +31 (0)30-2383749 nlinfo@phenomenex.com C18 00A-4435-E0 00D-4435-E0 Mexico t: 001-800-844-5226 C18 00A-4435-E0 00D-4435-E0 00D-4435-E0 t: 001-300-238-7768 tecnicomx@phenomenex.com * SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282 t: 030-2418700 Ordering Information Impact™ Precipitation Products	t: +31 (0)30-2418700 Phases 30 x 4.6 50 x 4.6 100 x 4.6 150 x 4.6 f: +31 (0)30-2383749 nlinfo@phenomenex.com C18 00B-4435-E0 00D-4435-E0 00F-4435-E0 Mexico t: 001-800-844-5226 * SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282 * SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282 The Netherlands t: 030-2418700 Ordering Information f: 030-2383749 Impact [™] Precipitation Products	t: +31 (0)30-2418700 Phases 30 x 4.6 50 x 4.6 100 x 4.6 150 x 4.6 250 x 4.6 f: +31 (0)30-2383749 nlinfo@phenomenex.com Mexico C18 00A-4435-E0 00B-4435-E0 00F-4435-E0 00G-4435-E0 t: 001-800-844-5226 * SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282 t: 030-2418700 * SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282 t: 030-2418700 Ordering Information f: 030-2383749 Impact [™] Precipitation Products

Canada

- t: (800) 543-3681
- f: (310) 328-7768
- info@phenomenex.com

Denmark

- t: 4824 8048
- f: +45 4810 6265 nordicinfo@phenomenex.com

beinfo@phenomenex.com

- Finland
- t: 09 4789 0063 f: +45 4810 6265
- nordicinfo@phenomenex.com

France

t: 01 30 09 21 10 f: 01 30 09 21 11 franceinfo@phenomenex.com

Germany

- t: 06021-58830-0
- f: 06021-58830-11 anfrage@phenomenex.com
 - India
- t: 040-3012 2400 f: 040-3012 2411 indiainfo@phenomenex.com

Ireland

- t: 01 247 5405
- f: +44 1625-501796 eireinfo@phenomenex.com

Italy

- t: 051 6327511 f: 051 6327555
- italiainfo@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

* SecurityGuard Analytical cartridges require holder, Part No.: KJ0-4282	
Ordering Information	

Impact[™] Precipitation Products

Part No.	Description	Unit
CE0-7565	Impact Protein Precipitation, Square Well, Filter Plate, 2 mL	2/box

Accessories

P

(

	-	
Collection	Plates (deep well, polypropylene)	
AH0-7192	Strata [®] 96-Well Collection Plate 350 µL/well	50/pk
AH0-7193	Strata 96-Well Collection Plate 1 mL/well	50/pk
AH0-7194	Strata 96-Well Collection Plate 2 mL/well	50/pk
AH0-8635	Strata 96-Well Collection Plate, 2 mL Square/Round-Conical	50/pk
	AH0-7192 AH0-7193 AH0-7194	AH0-7193 Strata 96-Well Collection Plate 1 mL/well AH0-7194 Strata 96-Well Collection Plate 2 mL/well



If Phenomenex products in this technical note do not provide at least an equivalent separation as compared to other products of the same phase and dimensions, return the product with comparative data within 45 days for a FULL REFUND.

Terms and Conditions Subject to Phenomenex Standard Terms and Conditions which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Gemini and Strata are registered trademarks of Phenomenex. Impact, Verex, and Security-Guard are trademarks of Phenomenex. API 4000 is a trademark of AB SCIEX Pte, Ltd. Gemini is patented by Phenomenex. U.S. Patent No. 7,563,367

SecurityGuard is patented by Phenomenex, U.S. Patent No. 6.162.362

© 2013 Phenomenex, Inc. All rights reserved.

() phenomenex