

The Evaluation of Multi-Pesticide Screening Methods by GC/MS

Sky Countryman, Jim Archer, Kory Kelly, and Doug Silva
Phenomenex, Inc., 411 Madrid Ave., Torrance, CA 90501 USA

Abstract

The current work demonstrates the use of two new and unique phases, which have been optimized for the analysis of all classes of pesticides. The phase chemistry improves separation and peak shape for the more polar pesticide compounds when compared to standard 5 % phenyl columns. Selectivity data is compared between a 5ms type phase and the two new columns.

Multi-pesticide residue screening has been evaluated for over 250 different pesticides commonly analyzed from fruits and vegetables (not all data presented here). The unique selectivity offered by the two phases improves resolution

for multi-component analytes providing a more unique elution pattern, which can be used to identify closely eluting analytes.

Since the phases have orthogonal selectivity, they are also a good choice for dual column methods. Some data is presented for EPA specified testing procedures.

Introduction

Pesticides are widely used by farmers to control pests, weeds, and molds that would otherwise decrease crop production. While this has significantly increased worldwide food productions, these same pesticides pose significant health and environmental risks. The restrictions for specific pesticides differ from one country to the next. As world trade increases, the potential threat to other countries' populations increases. For this reason, pesticides are the subjects of increasing regulation.

Since many different types of pesticides can be used on the same food product, Multi-Residue screening approaches are used to look for multiple classes of pesticide compounds at one time. Considering that there are more than 500 registered pesticides, no single analysis technique is capable of screening for all possible contaminants. However, gas chromatography (GC) is still the most commonly used method for the majority of the pesticide classes. While analyte specific detectors such as ECD or NPD may be used for screening, Mass Spectrometer (MS) detection must be employed to provide positive confirmation. Zebron MultiResidue™ (MR) columns were specially designed for pesticides analysis. The columns were developed using two new

stationary phases that are unlike any other commercially available columns. The phases were designed to provide orthogonal selectivity to provide maximum resolving power in complex samples. Zebron MultiResidue™ columns provide low bleed on ECD and NPD detectors and both columns are MS certified, so they can also be used with GC/MS for multi-residue pesticide methods.

GC/MS screening of multi-residue pesticide standards was evaluated using the new Zebron MultiResidue™ columns and compared with the results obtained using a standard 5ms type column. Dual column approaches were also evaluated using a chlorinated pesticides sample following EPA Method 8081A.

Figure 1: Multi-Pesticide Residue Analysis using Zebron MultiResidue™-1





Figure 2: Multi-Pesticide Residue Analysis using Zebron MultiResidue™-2



Peaks for Figure 1 and Figure 2:

Peak No. 1	Sample Analyte	Peak No. 1	Sample Analyte	Peak No. 1	Sample Analyte	Peak No. 1	Sample Analyte	Peak No. 1	Sample Analyte	Peak No. 1	Sample Analyte
1	Dichlorvos	18	Demeton	35	Phenothrin	62	Methyl parathion	89	MSK-CCA	101	Ethion
2	EPTC	19	Dichlorprop (methyl ester)	36	Atrazine	63	Phosphamidon isomer	90	Pendimethalin	102	Tricyclozole
3	Volinate	20	Trifluralin	37	Monosulfophos	64	Seckumion	91	Diphenamid	103	Permethrin
4	3,5-Dichlorobenzoic acid (methyl ester)	21	Thionex	38	Atrazine isomer	65	Dichlorfenthion	92	Carbendathion	104	Carbendathion
5	Veratrate	22	Cyfluthrin	39	Silver (methyl ester)	66	2,4-DB (methyl ester)	93	Butachlor	105	Famfur
6	Permethrin	23	Benflufen	40	Terbuthin	67	Terbuthin	94	Carbendathion	106	Nonfluorazone
7	Mevinphos	24	Propazine	41	Propazine	68	Dicofol (methyl ester)	95	Stirolos	107	Heptachlor
8	4-Methylphenol (methyl ester)	25	Chlorpyrifos	42	Alachlor	69	Alachlor	96	Triclorfon	108	EPN
9	Mevinphos isomer	26	Chlorpyrifos	43	Propanil	70	Propanil	97	Tridimenol	109	Phosmet
10	Trichlorfon	27	2,4-D (methyl ester)	44	Atrazine	71	Phosphamidon	98	Mesitox Goxide	110	Leptophos
11	Dicamba (methyl ester)	28	Sulfotop	45	Sinazine	72	Phosphamidon (methyl ester)	99	Neoprene	111	Azinphos-Methyl
12	MCPP (methyl ester)	29	Naled	46	Terbuthylazine	73	Prometryn	100	Fenarimol	112	Fenarimol
13	Molinate	30	Phorate	47	Disulfoton	74	Azinphos	101	Carbofent	113	Azinphos-ethyl
14	Tebuuthuron	31	Dicofol	48	Benfonos	75	Benfonos (methyl ester)	102	Carbofent	114	Coumaphos
15	MCPA (methyl ester)	32	Dimethoate	49	Dimethoate	76	Mesitox	103	Carbofent		
16	DEET	33	Profluralin	50	2,4,5-T Methyl ester	77	Sinetryn				
17	Terbuthyl pyrophosphate (methyl ester)	34	Demeton isomer	51	Disulfoton	78	Methidathion				

Figure 3: Nitrogen & Phosphorous Pesticides (NPM-102)




Figure 4: Nitrogen & Phosphorous Pesticides (NPM-106)




Figure 5: Nitrogen Pesticides (NPM-633)




Figure 6: Organophosphorous Pesticides (SPM-844)





Figure 7: Dual Column Organochlorine Pesticide Testing Under EPA Method 8081A



Results & Discussion

Analysis of 114 different pesticides was done using both Zebron MultiResidue-1 and MultiResidue-2 column (Figures 1 & 2). The Zebron MultiResidue-2 columns showed additional peaks for certain pesticides compared to the Zebron MultiResidue-1 column, however both columns provide increased resolution compared to a standard 5ms type phase.

Resolution of specific compounds was compared with the Zebron ZB-5ms phase to demonstrate the increased separation power offered by these new columns. In Figure 3, the elution order for atrazine and propazine is reversed on both the MultiResidue-1 and the MultiResidue-2 columns compared to the ZB-5ms. The MultiResidue-2 demonstrates much longer relative retention for these two compounds and is able to provide baseline resolution. In Figure 4, the MultiResidue-1 provides increased resolution of trifluralin and benflufen as well as several of the later eluting compounds. The polarity of the MultiResidue-2 caused the trifluralin to partially co-elute with propachlor, however resolution is substantially increased for the last four compounds.

Asymmetrical peak shape such as fronting is commonly observed for polar pesticides on non-polar phases because the pesticides are not soluble in the phase. This can lead to poor sensitivity and/or co-elutions with other closely eluting pesticides. Figure 5 shows overloading of terbutryn and bromacil on the ZB-5ms phase due to the non-polar character of the phase. The phase chemistry of the MultiResidue-1 and MultiResidue-2 was specially designed to match the polarity of pesticide compounds and provide better overall peak shape.

Many pesticides contain a mixture of isomers. Resolving the isomers can provide additional confirmation for samples with interferences in the

chromatographic region of the target analyte. In Figure 6, phosphamidon shows an isomer eluting just after terbutryn on the ZB-5ms. If high levels of terbutryn were present, it is unlikely resolution between the two compounds would be maintained. The MultiResidue-1 and the MultiResidue-2 columns significantly increase resolution of the phosphamidon isomer ensuring accurate identification.

The US EPA regulates the testing of 20 specific chlorinated pesticides under the official Method 8081A. The method specifies an Electron Capture Detector (ECD), which is extremely sensitive for chlorinated compounds, however, it does not provide any confirmatory information about the peak. To reduce the occurrence of misidentifications, the method requires the use of two GC columns of dissimilar selectivity in a parallel configuration. The EPA considers an analyte's presence confirmed if it has a peak at the pre-determined retention time on both columns.

The orthogonal selectivity provided by the Zebron MultiResidue-1 and the MultiResidue-2 allows for baseline resolution of all 20 chlorinated pesticides, surrogates, and internal standard in 10 minutes (Figure 7). Two elution order changes are observed between the phases, demonstrating that the phases are different enough to provide accurate confirmation.

Conclusion

As world trade increases, the potential threat to other countries' populations due to contaminated food products increases. Recent deaths caused by food exported from countries like China, emphasize the need for comprehensive testing procedures. The Zebron MultiResidue-1 and the MultiResidue-2 present a comprehensive solution for Multi-Pesticide residue testing by GC/MS and other hyphenated techniques. The columns provide additional confirmation, potential separation from matrix interferences, greater resolution of isomer peaks, and improved chromatography of more polar analytes.

Zebron and MultiResidue are trademarks of Phenomenex, Inc. © 2008 Phenomenex, Inc. All rights reserved.