

Pssst!
Waters® XBridge®
and **BEH™** user,
now you can...

Triumph

with

- Improved Peak Shape for Bases
- Shorter Run Times
- Core-Shell Performance from pH 1-12

Kinetex EVO C18



Core-Shell UHPLC/HPLC Columns



The Past

Traditional fully porous hybrid




- More band broadening
- Lower sensitivities
- Longer run times



The Present


Kinetex® Core-Shell EVO C18

- Sharper peaks
- Higher efficiencies
- Increased sensitivity
- Shorter run times



Introducing Kinetex EVO C18

The Kinetex core-shell EVO C18 provides much greater LC performance than fully porous products, giving you the freedom to **get improved results, cost savings, and overall higher productivity**. The patented surface technology of the Kinetex EVO C18 also provides you with the ruggedness to take on 100% aqueous conditions and reversed phase methods from pH 1-12.



Kinetex EVO C18 uses a patented organo-silica grafting process which incorporates uniform stabilizing ethane cross-linking to provide resistance to high pH attack while maintaining mechanical strength of the core-shell particle.

Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

© 2015 Phenomenex, Inc. All rights reserved.

How Does Kinetex EVO C18 Perform?



“With our vast range of clients, our quality control lab receives all kind of sample to analyse. This column is often used to analyse a large range of samples. The low cost combined with the great range of pH in which this column can be used saves us a lot of trouble and gives us a great price/utility ratio.”

Pierre-Alexandre Paquette
Solartech Laboratories, Inc.

“This column works better than any column I have used for method development.”

Jenn Jazikoff
Dynamic Pharmaceuticals

Upgrade from XBridge and BEH C18 columns:




XBridge® and BEH™ C18		Kinetex EVO C18	Page(s)
3.5 µm	vs	5 µm	4
5 µm	vs	5 µm	4-5
1.7 µm	vs	2.6 µm	6-7
1.7 µm	vs	1.7 µm	8-9
5 µm OBD™	vs	5 µm AXIA™	10
Ordering Information			11


Upgrade Your HPLC Methods for Greater Sensitivity

With the combination of 100% aqueous stability, pH ruggedness from 1-12, and the core-shell performance advantage, you can easily replace old hybrid silica columns and gain immediate method improvements without increasing backpressure.


Kinetex® 5µm EVO C18 150 x 4.6 mm



Waters® XBridge® 3.5µm C18 150 x 4.6 mm



Waters XBridge 5µm C18 150 x 4.6 mm



Conditions for all columns:

Column: Kinetex 5µm EVO C18
XBridge 5µm C18
XBridge 3.5µm C18

Dimensions: 150 x 4.6 mm

Mobile Phase: A: 0.1% Formic acid in Water
B: 0.1% Formic acid in Acetonitrile

Gradient: 5% to 95% B over 10 minutes

Flow Rate: 1.5 mL/min

Temperature: Ambient

Detection: UV @ 254 nm

Sample:


1. Uracil
2. Pindolol
3. Chlorpheniramine
4. Nortriptyline
5. 3-Methyl-4-Nitrobenzoic acid
6. 5-Methyl Salicyl Aldehyde
7. Hexaphenone

Comparative separations may not be representative of all applications.




Advance pH 1-12 Methods Ahead with Better Peak Shape

Easy improvements in resolution, productivity, and sensitivity are just a column away. By upgrading your methods from traditional fully porous hybrid materials to Kinetex 5 μ m EVO C18 Core-Shell Technology you can expect excellent performance and pH durability at low HPLC backpressures!

Kinetex 5 μ m EVO C18 150 x 4.6 mm



Waters® XBridge® 5 μ m C18 150 x 4.6 mm



Conditions for all columns:

Column: Kinetex 5 μ m EVO C18
XBridge 5 μ m C18

Dimensions: 150 x 4.6 mm

Mobile Phase: A: 20 mM Sodium Phosphate (pH 7.0 or pH 10.2)
B: Methanol/Acetonitrile (50:50)

Isocratic: A/B (30:70)

Flow Rate: 1.25 mL/min

Temperature: 30 °C

Detection: UV @ 254 nm

Sample: 1. Protriptyline
2. Nortriptyline
3. Imipramine
4. Amitriptyline
5. Clomipramine

Wow! Better peak shape and better performance in less time!







Replace Fully Porous Sub-2 μ m Columns for Greater Productivity


Utilize Kinetex[®] 2.6 μ m core-shell columns to take full advantage of the flow rates that are available on your UHPLC system. Since Kinetex 2.6 μ m core-shell columns display similar or better performance at a lower pressure than fully porous sub-2 μ m columns, you can increase your flow rate with the Kinetex 2.6 μ m column and achieve some very useful time savings and productivity gains, while holding onto the necessary performance parameters.

Waters[®] ACQUITY[®] BEH[™] 1.7 μ m C18



35% Faster Analysis

Kinetex 2.6 μ m EVO C18




Conditions for all columns unless otherwise noted:

Columns: Kinetex 2.6 μ m EVO C18	Flow Rate: 0.8 mL/min (Kinetex)
ACQUITY BEH 1.7 μ m C18	0.4 mL/min (BEH)
Dimensions: 50 x 2.1 mm	Temperature: 30 °C
Mobile Phase: A: Water with 0.1% Formic acid	Detection: MS/MS
B: Acetonitrile with 0.1% Formic acid	Sample: Proprietary Drug Formulation
Gradient:	
Time (min) %B	
0 5	
5 95	


Comparative separations may not be representative of all applications.

Replace Fully Porous Sub-2 μ m Columns for Increased Resolution

Start Run a Kinetex 2.6 μ m column in the same dimension as your sub-2 μ m fully porous column.



What's Next? Utilize a longer Kinetex column to further increase resolution



Conditions for all columns unless otherwise noted:


Column: Kinetex 2.6 μ m EVO C18 ACQUITY BEH 1.7 μ m C18	Flow Rate: 0.4 mL/min
Dimensions: 50 x 2.1 mm	Temperature: 22 °C
Mobile Phase: A: 20 mM Potassium phosphate (pH 2.5) B: Acetonitrile	Detection: UV @ 210 nm
Gradient:	Sample: 1. Green Coffee Extract (Chlorogenic Acid and Antioxidants)
Time (min)	% B
0	5
5	70
6	70
6.5	5
10	5

Comparative separations may not be representative of all applications.

Improve Your UHPLC Analyses with Faster Results and Better Sensitivity

By upgrading from fully porous hybrid UHPLC columns to the Kinetex® 1.7 µm EVO C18, you can keep the pH stability but gain much greater performance. With the core-shell advantage at your disposal your new methods will be faster and more sensitive, leaving room for greater productivity.

Berberine Standard



57% Increase in Peak Height Response

After enjoying the greater sensitivity, resolution, and peak shape given by the Kinetex EVO C18 for this isolated standard, the same performance benefits can easily be realized for the more complex extracted sample.


Conditions for both columns:

Column: Kinetex 1.7 µm EVO C18
ACQUITY BEH 1.7 µm C18
Dimensions: 50 x 2.1 mm
Mobile Phase: A: 0.1 % Formic Acid in Water
B: 0.1 % Formic Acid in Methanol

Gradient:	Time (min)	% B
	0	2
	5	95

Flow Rate: 0.5 mL/min
Temperature: 30°C
Detection: UV @ 254 nm
Sample: 1. Berberine

Goldenseal Extract (Containing Berberine)




Comparative separations may not be representative of all applications.

More Performance, Excelling Quality

Here's an illustration of just how powerful a tool the Kinetex® 1.7 µm EVO C18 is. In comparison with multiple batches of a common hybrid sub-2 µm column line, Kinetex EVO C18 columns show consistently greater performance and better peak shape for basic compounds. This level of quality is what we strive to give you every day with every batch, every column, and even every phone call.

Kinetex 1.7 µm EVO C18



Conditions for both columns:

Column: Kinetex 1.7 µm EVO C18
ACQUITY BEH 1.7 µm C18

Dimensions: 150 x 2.1 mm

Mobile Phase: A: 0.1 % Formic Acid in Water
B: 0.1 % Formic Acid in Acetonitrile

Gradient:	Time (min)	% B
	0	5
	5	95
	6	95
	6.5	5
	10	5


Flow Rate: 0.5 mL/min

Temperature: 22 °C


Detection: UV @ 254 nm

Sample: 1. 4-Hydroxybenzoic Acid
2. Labetolol
3. Propranolol
4. Protriptyline
5. Ibuprofen

Average Peak Capacity_{A&B}: 262




Waters® ACQUITY® BEH™ 1.7 µm C18



Wow! Look at my improved peak shape for bases!

Average Peak Capacity_{A&B}: 168




Comparative separations may not be representative of all applications.


Enhance Purification and Isolation Levels

Kinetex® EVO C18 will easily pave the way to enhanced purification and isolation levels with its incredible pH durability and core-shell performance gains. Get ready for a future of high yields and chromatographic improvements over traditional fully porous and hybrid materials.

Kinetex 5µm EVO C18
150 x 21.2mm AXIA™



Waters® XBridge® 5µm C18
150 x 19mm Prep OBD™



Conditions for all columns:

- Column:** Kinetex 5µm EVO C18 Axia packed
XBridge 5µm C18 Prep OBD
- Dimensions:** 150 x 21.2mm (Kinetex)
150 x 19mm (XBridge)
- Mobile Phase:** A: 0.1 % TFA in Water
B: 0.1 % TFA in Acetonitrile
- Gradient:** 5-95 % B in 4 min.
Hold for 1 min.
Re-equilibrate for 2 min.
- Flow Rate:** 25 mL/min
- Temperature:** Ambient
- Detection:** UV @ 254 nm
- Sample:** Amitriptyline


- Sample Load:**
- 2 mg
 - 4 mg
 - 10 mg
 - 20 mg
 - 40 mg
 - 80 mg

Learn more online about how Kinetex EVO C18 can enhance your HPLC methods!

[www.phenomenex.com/
KinetexEVO](http://www.phenomenex.com/KinetexEVO)



[www.phenomenex.com/
KinetexApps](http://www.phenomenex.com/KinetexApps)



Comparative separations may not be representative of all applications.

Ordering Information

5 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges [‡]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4633-AN	00B-4633-AN	00D-4633-AN	00F-4633-AN	AJ0-9298 for 2.1 mm ID

5 µm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	00B-4633-Y0	00D-4633-Y0	00F-4633-Y0	AJ0-9297 for 3.0 mm ID

5 µm Analytical Columns (mm)					SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 4.6	100 x 4.6	150 x 4.6	250 x 4.6	3/pk
EVO C18	00B-4633-E0	00D-4633-E0	00F-4633-E0	00G-4633-E0	AJ0-9296 for 4.6 mm ID

5 µm Axia™ Packed Preparative Columns (mm)					SecurityGuard PREP Cartridges*
Phases	50 x 21.2	100 x 21.2	150 x 21.2	250 x 21.2	15 x 21.2
EVO C18	00B-4633-P0-AX	00D-4633-P0-AX	00F-4633-P0-AX	00G-4633-P0-AX	AJ0-9304 for 21.2 mm ID

5 µm Axia Packed Preparative Columns (mm)				SecurityGuard PREP Cartridges**
Phases	50 x 30	100 x 30	250 x 30	15 x 30
EVO C18	00B-4633-U0-AX	00D-4633-U0-AX	00G-4633-U0-AX	AJ0-9305 for 30 mm ID

2.6 µm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges [‡]
Phases	30 x 2.1	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00A-4725-AN	00B-4725-AN	00D-4725-AN	00F-4725-AN	AJ0-9298 for 2.1 mm ID

2.6 µm MidBore™ Columns (mm)				SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 3.0	100 x 3.0	150 x 3.0	3/pk
EVO C18	00B-4725-Y0	00D-4725-Y0	00F-4725-Y0	AJ0-9297 for 3.0 mm ID

2.6 µm Analytical Columns (mm)				SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 4.6	100 x 4.6	150 x 4.6	3/pk
EVO C18	00B-4725-E0	00D-4725-E0	00F-4725-E0	AJ0-9296 for 4.6 mm ID

1.7 µm Minibore Columns (mm)				SecurityGuard ULTRA Cartridges [‡]
Phases	50 x 2.1	100 x 2.1	150 x 2.1	3/pk
EVO C18	00B-4726-AN	00D-4726-AN	00F-4726-AN	AJ0-9298 for 2.1 mm ID

[‡]SecurityGuard ULTRA Cartridges require holder, Part No.: AJ0-9000

*PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8223

**PREP SecurityGuard Cartridges require holder, Part No.: AJ0-8277



If you are not completely satisfied with Kinetex core-shell columns, send in your comparative data to a similar product with the Kinetex column within 45 days for a FULL REFUND.

Triumph with

- Improved Peak Shape for Bases
- Shorter Run Times
- Core-Shell Performance from pH 1-12

Kinetex EVO C18 Core-Shell UHPLC/HPLC Columns

Australia

t: +61 (0)2-9428-6444
f: +61 (0)2-9428-6445
auinfo@phenomenex.com

Austria

t: +43 (0)1-319-1301
f: +43 (0)1-319-1300
anfrage@phenomenex.com

Belgium

t: +32 (0)2 503 4015 (French)
t: +32 (0)2 511 8666 (Dutch)
f: +31 (0)30-2383749
beinfo@phenomenex.com

Canada

t: +1 (800) 543-3681
f: +1 (310) 328-7768
info@phenomenex.com

China

t: +86 (0)20 2282-6668
f: +86 (0)20 2809-8130
chinainfo@phenomenex.com

Denmark

t: +45 4824 8048
f: +45 4810 6265
nordicinfo@phenomenex.com

Finland

t: +358 (0)9 4789 0063
f: +45 4810 6265
nordicinfo@phenomenex.com

France

t: +33 (0)1 30 09 21 10
f: +33 (0)1 30 09 21 11
franceinfo@phenomenex.com

Germany

t: +49 (0)6021-58830-0
f: +49 (0)6021-58830-11
anfrage@phenomenex.com

India

t: +91 (0)40-3012 2400
f: +91 (0)40-3012 2411
indiainfo@phenomenex.com

Ireland

t: +353 (0)1 247 5405
f: +44 1625-501796
eireinfo@phenomenex.com

Italy

t: +39 051 6327511
f: +39 051 6327555
italiainfo@phenomenex.com

Luxembourg

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

Mexico

t: 001-800-844-5226
f: 001-310-328-7768
tecnicomx@phenomenex.com

The Netherlands

t: +31 (0)30-2418700
f: +31 (0)30-2383749
nlinfo@phenomenex.com

New Zealand

t: +64 (0)9-4780951
f: +64 (0)9-4780952
nzinfo@phenomenex.com

Norway

t: +47 810 02 005
f: +45 4810 6265
nordicinfo@phenomenex.com

Puerto Rico

t: +1 (800) 541-HPLC
f: +1 (310) 328-7768
info@phenomenex.com

Spain

t: +34 91-413-8613
f: +34 91-413-2290
espinfo@phenomenex.com

Sweden

t: +46 (0)8 611 6950
f: +45 4810 6265
nordicinfo@phenomenex.com

United Kingdom

t: +44 (0)1625-501367
f: +44 (0)1625-501796
ukinfo@phenomenex.com

USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com

All other countries Corporate Office USA

t: +1 (310) 212-0555
f: +1 (310) 328-7768
info@phenomenex.com



www.phenomenex.com

Phenomenex products are available worldwide. For the distributor in your country, contact Phenomenex USA, International Department at international@phenomenex.com

Terms and Conditions

Subject to Phenomenex Standard Terms & Conditions, which may be viewed at www.phenomenex.com/TermsAndConditions.

Trademarks

Phenomenex and Kinetex are registered trademarks of Phenomenex. Axia, MidBore, and SecurityGuard are trademarks of Phenomenex. Waters, ACQUITY, UPLC, and XBridge are registered trademarks, BEH Technology and OBD are trademarks of Waters Corporation. SelectScience is a registered trademark of SelectScience Ltd.

Disclaimer

Phenomenex is not affiliated with Waters Corporation. Comparative separations may not be representative of all applications.

Axia column and packing technology is patented by Phenomenex. U.S. Patent No. 7, 674, 383 Kinetex EVO is patented by Phenomenex. U.S. Patent Nos. 7,563,367 and 8,658,038 and foreign counterparts.

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

© 2015 Phenomenex, Inc. All rights reserved.